

TECHNICAL INFORMATION

Catalog Number: 150025

Fibronectin

Molecular Weight: 220,000

Physical Description: White lyophilized powder

Source: Bovine Plasma


Description: Fibronectins are large glycoproteins found in cells and extracellular fluids.¹ They are involved in a wide range of biological activities including cell adhesion to extracellular matrices^{2,3}, blood coagulation^{4,5} and wound healing.⁶ In tissue culture, fibronectins promote cell adhesion and spreading on synthetic surfaces.^{7, 8, 9} Certain cell types show morphological changes relating to loss of fibronectin from the cell surface^{10,11,12,13} or depletion of fibronectin from the culture medium.¹⁴ Fibronectins are of two, closely related, major types¹⁵: cellular fibronectin, otherwise known as Large External Transformation Sensitive (LETS) protein¹⁶ and plasma fibronectin, also called Cold-insoluble Globulin (CLG).¹⁷

Fibronectin is present in plasma at a concentration of 250-300 ug/mL. When plasma is stored at +4°C, a cryoprecipitate is formed containing large amounts of fibronectin. This property gave rise to the early name of "cold-insoluble globulin".¹⁷ Plasma Fibronectin exists in multiple molecular forms, 90% of which is a polypeptide dimer of 220,000 daltons.¹⁸ Sodium dodecyl sulfate/polyacrylamide (SDS-PAGE) gel electrophoresis of reduced plasma fibronectin resolves the dimer as two closely spaced bands. Under non-reducing conditions, an additional band is observed corresponding to a molecular weight of 450,000 daltons.

Properties of Fibronectin¹⁹

pI	5.5-6.3
Carbohydrate	5-9% asparagine-linked complex oligosaccharides
Sulphydryl Groups	One or two in the C-terminal 30%
Disulphide Bonds	~20 per sub unit. Inter sub-unit bond(s) are very near C-terminal. N terminal is very rich in interchain disulphides.
Tertiary Structure	No a helix, probably some b structure
Secondary Structure	Asymmetric and elongated with globular domains

Functional Domains of Fibronectin:¹⁹



The regions depicted by thin lines are extremely sensitive to proteolysis. Other regions (I-VI) are compact and globular and contain specific binding sites for other molecules.

Biological Activities of Fibronectin:

The biological activities of fibronectin are largely associated with its property of being able to interact with many different substrates and molecules. Some of these properties are listed in the table below.

Activity	Reference number
Cell-Cell aggregation	20
Live cells	21
Erythrocytes	
Cell-substratum adhesiveness	22
Collagen substratum	7,8
Plastic substratum	23
Cell spreading	

Reversion of Transformed Phenotype	24
Cell alignment and overlapping	13,14
Fibroblastic morphology	11,12
Cell surface morphology	25
Microfilament bundle organization	
Increased Cell Motility	26
Phagocytosis	27
Embryonic Differentiation	28,29

Uses of Fibronectin:

Fibronectin has wide applications in cellular research into proliferation, differentiation, development and morphology. Precoating of glass, plastic and micro-carrier beads aids cell attachment, particularly those cells which are deficient in cell surface fibronectin. It has been reported that proteolytic digests of fibronectin are mitogenic for CH23 fibroblasts.³⁰

In summary, fibronectin will:

- Bind to commonly used cell-culture surfaces.
- Increase plating efficiency of fastidious cells.
- Increase cloning efficiency.
- Reduce variation in plating efficiency found in different serum batches.

Suggested Procedure for Coating Cell Cultureware

- Determine the amount of fibronectin needed to coat culture vessels by multiplying the total surface area (cm²) by the desired concentration (mg/mL) of fibronectin. Recommended amount is 2-10 mg/cm².
- Reconstitute fibronectin at 1 mg/mL and dilute to the desired concentration with physiological buffer. Wet the surface of each culture vessel to be coated with a minimum amount of sterile balanced salt solution (serum and protein free) required to cover the entire area.
- Introduce the proper CO₂ atmosphere, if required.
- Add the calculated amount of fibronectin to each culture vessel.
- Allow fibronectin to adsorb to the surface of the vessel for 5-20 minutes.
- Remove residual balanced salt solution before proceeding with standard cell culture procedures.

Solubility: Soluble in water

References:

- Yamada, K.M., and Olden, K., "Fibronectins - adhesive glycoproteins of cell surface and blood." *Nature*, **v. 275**, 179-184 (1978).
- Grinnell, F., "Cellular adhesiveness and extracellular substrata." *Int. Rev. Cytol.*, **v. 53**, 65-144 (1978).
- Grinnell, F., Feld, M., and Minter, D., "Fibroblast adhesion to fibronectin and fibrin substrata: requirement for cold-insoluble globulin (plasma fibronectin)." *Cell*, **v. 19**, 517-525 (1980).
- Arneson, M.A., Hammerschmidt, D.E., Furcht, L.T., and King, R.A., "A new form of Ehlers-Danlos Syndrome. Fibronectin corrects defective platelet function." *JAMA*, **v. 244**, 144-147 (1980).
- Mosher, D.F., Schad, P.E., and Vann, J.W., "Cross-linking of collagen and fibronectin by factor XIIIa: localisation of participating glutamyl residues to a tryptic fragment of fibronectin." *J. Biol. Chem.*, **v. 255**, 8848-8858 (1980).
- Grinnell, F., Billingham, R.E., and Burgess, L., "Distribution of fibronectin during wound healing in vivo." *J. Invest. Dermatol.*, **v. 75**, 181-189 (1981).
- Norton, E.K., and Izzard, C.S., "Fibronectin promotes formation of the close cell-to-substrate contact in cultured cells." *Exp. Cell Res.*, **v. 139**, 463-467 (1982).
- Hook, M., Rubin, K., Oldberg, A., Obrink, B., and Vaheri, A., "Cold-insoluble globulin mediates the adhesion of rat liver cell to plastic petri dishes." *Biochim. Biophys. Res. Commun.*, **v. 79**, 726-731 (1977).
- Schwarz, M.A., and Juliano, R.L., "Interaction of fibronectin-coated beads with CHO cells." *Exp. Cell Res.*, **v. 152**, 302-312 (1984).
- Hynes, R.O., "Cell surface proteins and malignant transformation." *Biochim. Biophys. Acta.*, **v. 458**, 73-107 (1976).
- Hynes, R.O., Ali, I.U., Destree, A.T., Mautner, V.M., Perkins, M.E., Senger, D.R., Wagner, D.D., and Smith, K.K., "A large glycoprotein lost from the surfaces of transformed cells." *Ann. N.Y. Acad. Sci.*, **v. 312**, 317-342 (1978).
- Hayman, E.G., Engvall, E., and Ruoslahti, E., "Concomitant loss of cell surface fibronectin and laminin from transformed rat kidney cells." *J. Cell Biol.*, **v. 88**, 353-357 (1981).
- Yamada, K.M., Olden, K., and Pastan, I., "Transformation - sensitive cell surface protein: isolation, characterization, and role in cellular morphology and adhesion." *Ann. N.Y. Acad. Sci.*, **v. 312**, 256-277 (1978).
- Ali, I.U., Mautner, V.M., Lanza, R.P. and Hynes, R.O., "Restoration of normal morphology, adhesion and cytoskeleton in transformed cells by addition of a transformation-sensitive surface protein." *Cell*, **v. 11**, 115-126 (1977).
- Yamada, K.M., and Kennedy, D.W., "Fibroblast cellular and plasma fibronectins are similar but not identical." *J. Cell Biol.*, **v. 80**, 492-498 (1979).
- Hynes, R.O., and Bye, J.M., "Density and cell cycle dependence of cell surface proteins in hamster fibroblasts." *Cell*, **v. 3**, 113-120 (1974).
- Morrison, P., Edsall, R., and Miller, S.G., "Preparation and properties of serum and plasma protein XIII. The separation of purified fibrinogen from fraction I of human plasma." *J. Am. Chem. Soc.*, **v. 70**, 3103-3108 (1948).
- Amrani, D.L., Homandberg, G.A., Tooney, N.M., Wolfenstein-Todel, C., and Mosesson, M.W., "Separation and analysis of the major forms of plasma fibronectin." *Biochim. Biophys. Acta.*, **v. 748**, 308-320 (1983).
- Hynes, R.O., and Yamada, K.M., "Fibronectins: Multifunctional modular glycoproteins." *J. Cell Biol.*, **v. 95**, 369-377 (1982).

- Perkins, M.E., Ji, T.H., and Hynes, R.O., "Cross linking of fibronectin to sulphated proteoglycans at the cell surface." *Cell*, **v. 16**, 941-952 (1979).
- Yamada, K.M., Yamada, S.S., and Pastan, I., "The major cell surface glycoprotein of chick embryo fibroblasts is an agglutinin." *Proc. Natl. Acad. Sci. (Wash.)*, **v. 72**, 3158-3162 (1975).
- Kleinman, H.K., Klebe, R.J., and Martin, G.R., "Role of collagenous matrices in the adhesion and growth of cells." *J. Cell Biol.*, **v. 88**, 473-485 (1981).
- Grinnell, F., and Minter, D., "Attachment and spreading of baby hamster kidney cells to collagen substrata: effects of cold-insoluble globulin." *Proc. Natl. Acad. Sci. (USA)*, **v. 75**, 4408-4412 (1977).
- Dodson, M.G., Gelder, F.B., Slota, J., and Lange, C., "Fibronectin and anchorage-independent and anchorage-dependent growth of benign and malignant cell lines." *Int. J. Cancer*, **v. 32**, 211-217 (1983).
- Keski-Oja, J., Sen, A., and Todaro, G.J., "Direct association of fibronectin and actin molecules in vitro." *J. Cell Biol.*, **v. 85**, 527-533 (1980).
- Gauss-Muller, V., Kleinman, H.K., Martin, G.R., and Schiffman, E., "Role of attachment and attractants in fibroblast chemotaxis." *J. Lab. Clin. Med.*, **v. 96**, 1071-1080 (1980).
- van de Water III, L., Schroeder, S., Grenshaw, E.B., and Hynes, R.O., "Phagocytosis of gelatin-latex particles by a murine macrophage line is dependent on fibronectin and heparin." *J. Cell Biol.*, **v. 90**, 32-39 (1981).
- Heasman, J., Hynes, R.O., Swan, A.P., Thomas, V., and Wylie, C.C., "Primordial germ cells of *Xenopus* embryos: the role of fibronectin in their adhesion during migration." *Cell*, **v. 27**, 437-447 (1981).
- Katow, H., Yamada, K.M., and Solursh, M., "Occurrence of fibronectin on the primary mesenchyme cell surface during migration in the sea urchin embryo." *Differentiation*, **v. 28**, 120-124 (1982).
- Humphries, M.J., and Ayad, S.R., "Stimulation of DNA synthesis by cathepsin-D digests of fibronectin." *Nature*, **v. 305**, 811-813 (1983).