



MP Biomedicals, LLC

29525 Fountain Parkway
Solon, Ohio 44139

Telephone: 440/337-1200
Toll Free: 800/854-0530
Fax: 440/337-1180
mailto: biotech@mpbio.com
web: <http://www.mpbio.com>

TECHNICAL INFORMATION

Catalog Number: 2030046, 2030026, 20300M2

CellVation Cryopreservation Medium

CellVation™: An alternative serum-free cryopreservation medium

To complete the serum-free product range, MP introduces CellVation™, a unique cryopreservation medium which does not contain either DMSO or serum. This cell freezing medium offers the following benefits:

- No presence of serum during freeze/thaw-cycles of cells adapted to serum-free growth
- Absence of cytotoxicity or morphological differentiation associated with DMSO
- Optimal recovery and faster cell attachment to substrate
- Ready to use

CellVation™ has a proven track record of successful freeze/thaw cycles of over 30 cell types, including primary and finite life cells as well as suspension and anchorage-dependent cell lines.

CellVation™ Cryopreservation Medium is

- Completely defined - No serum
- Ready-to-use solution
- Low cytotoxicity: No DMSO (Dimethyl sulfoxide)
- Proven track record of successful freezing of > 50 cell types: cell lines and primary cells
- Completes the MP serum-free product range: cells that have been adapted to serum free growth can then be stored in serum-free environment.

The following is a brief listing of the cell types that have been successfully frozen in MP's cryopreservation medium CellVation™.

Cell types used with CellVation™

Cell Type Description

A-375 Malignant Melanoma, Human
A549 Lung Carcinoma, Human
D283 Med Medulloblastoma, Human
Daoy Medulloblastoma, Human
DU 145 Carcinoma, Prostate, Metastasis to Brain, Human
HeLa Epitheloid Carcinoma, Cervix, Human
Hep-2 Epidermoid Carcinoma, Larynx, Human
LNCaP Metastatic Prostate Adenocarcinoma, Human
LS 174T Colon, Adenocarcinoma, Human
MCF7 Breast Adenocarcinoma, Pleural Effusion, Human
MCF 10 Mammary Gland, Human
MRC-5 Embryonal Lung, Diploid, Male, Human
NCI-H69 Small Cell Carcinoma, Lung, Human
PC-3 Prostate Adenocarcinoma, Human
Primary Bone Marrow and Peripheral Blood Stem Cells, Human
Primary Lymphocytes, Human
Raji Burkitt's Lymphoma, Human
SF Foreskin, Human
ST486 Burkitt's Lymphoma, Human
U-937 Histiocytic Lymphoma, Human
3T3 Embryonic Fibroblast, Mouse
B16BL6 Mouse, Melanoma
B16F10 Mouse, Melanoma
NSOMyeloma, Mouse
P31NSI/Ag4-1(NS-1) Non-secreting Myeloma, Mouse
Sp2/0-Ag14 Hybridoma, Non-secreting, mouse
Rat2 Embryo, Thymidine Kinase Mutant, Rat
CBAorta Endothelial, Bovine (primary)
MDBK Kidney, Bovine

CHO-K1 Ovary, Chinese Hamster
 hiDCK Kidney, Canine
 CRFK Kidney, Feline
 LLC-PK Kidney, Porcine
 STFetal Testis, Porcine (primary)
 Vero Kidney, African Green Monkey
 ZF4 Embryonic Fibroblastic Cells, Zebra Fish
 JB6C1 30-7b Epidermis Tumor, Mouse
 JB6C1 41-5a Epidermis Tumor, Mouse
 RT 101 Epidermis Tumor, Mouse
 RAW 264.7 Monocyte-Macrophage, Mouse
 CPA-47 Pulmonary Artery, Bovine
 CPAE Pulmonary Artery Endothelium, Bos Taurus
 PU5-1.8 (PU5-1R) Monocyte-Macrophage, Mouse
 IC-21 Macrophage, Mouse
 P388D Monocyte-Macrophage, Mouse
 BL-3B Lymphocytes, Leukemia Cell Suspension
 BTTurbinate Cells, Bovine, Bos Taurus
 CCD-18Co Colon Fibroblast, Normal, Human
 EBTr Embryonic Trachea, Bovine, Bos Taurus
 EJG Capillary Endothelium, Bovine
 HCT 116 Colon Carcinoma, Human
 LB9.K Kidney, Normal, Bovine
 6LBNL Lymph node, Ca., Bovine
 LS180 Colon, Adenocarcinoma, Human
 293 Transformed Primary Embryonal Kidney, Human

Typical Recovery Rate:

Cell Line	% Viability before freezing	% Viability after freezing with CellVation	% Viability after freezing with 10% DMSO, 10% FBS, RPMI 1640
SP 2/0 AG 14 Mouse Myeloma	100%	95%	90%
A 375 Human Malignant Melanoma	95%	85%	80%
U937 Human Histiocytic Lymphoma	90%	85%	60%

TCM A general purpose serum replacement 2 X 10 ml 2010026
 produces 1 L of medium
 2 X 100 ml 2010022
 produces 10 L of medium
 TCH Does not contain animal proteins, primarily 2 X 10 ml 220026
 for human cell cultures. Effective also on produces 1 L of medium
 other mammalian cells
 2 X 100 ml 2020022
 produces 10 L of medium
 CellVation™ cryopreservation medium without 20 ml 2030046
 serum and DMSO 60 ml 20300M2

Procedure for Use

CellVation™ is a cryopreservation medium which does not contain DMSO or serum. CellVation™ provides optimal recovery and is convenient to use. CellVation™ is packaged as a 1X concentrate and should not be diluted.

Recommended Freezing and Thawing Protocol

Freezing: Mix well before use. Do not dilute.

- Examine the culture for the absence of contamination, healthy growth, confluency, etc.
- If freezing adherent cells, remove using 0.25% trypsin for 1-3 minutes at 37°C. NOTE: Some cell lines grown in serum-free medium may be sensitive to 0.25% trypsin and therefore, may require less trypsin or the addition of a trypsin inhibitor.
- Perform a cell count to determine the total number of viable cells. Cell viability should be greater than 80% and cells should be in late log phase or pre-confluency growth phase.
- Centrifuge cells at 600-800 rpms for 10 minutes. Remove supernatant and save 3-5 ml for sterility testing (e.g. thioglycollate, brain heart infusion, etc.) and mycoplasma testing.
- Resuspend the cells gently in an appropriate volume of CellVation™ at a concentration of 1 x 10⁶ - 1 x 10⁷ cells/ml. Some cell types such as hybridomas and myelomas may require an increase in cell density.
- Dispense the cell suspension in 1-2 ml aliquots in plastic or glass ampoules.
- Seal ampoules and store at room temperature for 20 - 30 minutes.
- Place ampoules in an insulated container and store in a -70°C freezer for two hours. Remove insulation and transfer to the vapor phase of liquid nitrogen. If a -70°C freezer is not available, place an insulated container in the vapor phase of the liquid nitrogen. Remove the insulation after two hours and store vials in vapor phase for 24 hours before transferring to liquid phase. The suggested optimum cooling rate is 1°C per minute.

Recovery

- Remove vials from freezer and rapidly thaw in a 37°C waterbath.
- Wipe vials with 70% ethanol.
- Transfer cells to a culture flask and slowly add the appropriate volume of growth medium.
- Cell viability count should be performed at least 2 hours after recovery.

If desired, cells may be centrifuged in the following manner.

- Transfer cells to a 15 ml centrifuge tube and slowly add 2 - 3 mls of complete growth medium.
- Centrifuge at 400 - 600 rpms for approximately 5 minutes.
- Decant, resuspend the cells in the appropriate volume of growth medium, and transfer to culture flask.
- Cell viability count should be performed at least two hours after recovery.

IMPORTANT: Before terminating a culture, it is recommended to complete an entire freeze/thaw cycle to ensure sterility of the culture and cell viability.

Reference:

Pirro, J.P., et. al., *J. of Nuclear Medicine*, v. **35**, #9 (Sept. 1994)

Any questions concerning MP products can be addressed directly to our technical services department: 1-800-854-0530 or 1-330-447-1200.

FOR RESEARCH USE ONLY. NOT FOR HUMAN USE.