

4

Diagnostic Specificity (In-house study)

Specificity performance of individual markers

Method used	Reference assay: commercial ELISAs (n=400)#		
	Markers	Negative	Specificity
MULTISURE® Dengue Ab/ Ag Rapid Test	IgG	393	98.25%
	IgM	386	96.50%
	IgA	398	99.50%
	NS1	400	100.00%
	All markers	379	94.75%

Dengue negative samples are defined by Commercial dengue IgG, IgM and IgA ELISAs.

Specificity performance in various populations

Sample Category	Total Sample Size	Diagnostic Specificity of MULTISURE® Dengue Ab/Ag Rapid Test
Blood Donors / Healthy Donors	153	96.73% (148/153)
Hospitalized / Clinical	91	94.51% (86/91)
Cross-reactive	62	96.77% (60/62)
Interference	50	98.00% (49/50)
Pregnancy	44	81.82% (36/44)
Total	400	94.75% (379/400)

Specificity performance in potential cross-reactive specimens

Sample Profile	Total Sample Size	Performance of MULTISURE® Dengue Ab/Ag Rapid Test
Rheumatoid Factor (RF)	10	100.00% (10/10)
Hepatitis A	10	100.00% (10/10)
HBsAg	3	100.00% (03/03)
Hepatitis C	6	100.00% (06/06)
Measles IgG	9	100.00% (09/09)
Human T-lymphotropic Virus Type I (HTLV-I)	10	90.00% (09/10)
HIV 1/2	9	88.89% (08/09)
Hepatitis E	5	100.00% (05/05)
Total	62	96.77% (60/62)

Specificity performance in potential interference specimens

Sample Profile	Total Sample Size	Performance of MULTISURE® Dengue Ab/Ag Rapid Test
Icteric	8	100.00% (08/08)
Hemolyzed	8	100.00% (08/08)
Triglyceride	8	100.00% (08/08)
Lipemic	6	100.00% (06/06)
Total Protein	4	100.00% (04/04)
Total bilirubin	7	100.00% (07/07)
Antinuclear Antibody (ANA)	9	88.89% (08/09)
Total	50	98.00% (49/50)

External evaluation study

Sensitivity performance on different types of dengue infection

Dengue Type	MULTISURE® Dengue Ab/Ag Rapid Test (n=80)			
	Positive	Negative	Total	Sensitivity
DENV-1	17	3	20	85.0%
DENV-2	16	4	20	80.0%
DENV-3	20	0	20	100.0%
DENV-4	19	1	20	95.0%
Total	72	8	80	90.0%

All specimens were confirmed by PCR test.

Specificity performance on cross-reactivity test with other mosquito-borne diseases

Sample Profile	Total Sample size	Specificity
Chikungunya	5	100.0% (5/5)
Japanese Encephalitis	9	88.9% (8/9)
Total	14	92.9% (13/14)

Cross-reactivity with other flavivirus mediated and mosquito-borne diseases are known to be common but have not been tested.

Analytical Sensitivity

Analytical sensitivity for MP Diagnostics MULTISURE® Dengue Ab/Ag Rapid Test was found to be comparable to a leading commercial dengue rapid test for antibodies and NS1 antigen detection.

Precision

The inter-assay (between-run) and intra-assay (within-run, within-day and day-to-day) reproducibility of MP Diagnostics MULTISURE® Dengue Ab/Ag Rapid Test have been evaluated using a set of control panel members. All results obtained were consistently fall within the acceptance criteria, indicating the MP Diagnostics MULTISURE® Dengue Ab/Ag Rapid Test is robust, reproducible and consistent across three lots studied.

REFERENCES

1. Blacksell SD. Review Article: Commercial dengue rapid diagnostic test for point-of-care application: recent evaluations and future needs? J Biomed Biotechnol. 2012;2012:151967. doi: 10.1155/2012/151967. Epub 2012 May 10.
2. Weiss RC, Scott FW. Antibody-mediated enhancement of disease in feline infectious peritonitis: Comparative with dengue hemorrhagic fever. Comparative Immunology, Microbiology and Infectious Diseases 1981; Vol. 4, 2:175-189.
3. Guzman MG, Scott BH, Harvey A, Philippe B, Jeremy F, Duane JG, Elizabeth H, Axel K, Harold SM, Eric M, Michael BN, Jose LP, Cameron S, Sutee Y, Rosanna WP. Dengue: a continuing global threat. WHO/ TDR, 2010. Nature reviews: Microbiology. December 2010, S7-16.
4. Alvina CF, Camila MR, Cristiane CC, Célia LR, Lucy VB, Evaldo SA, Andréia MM, Karina IC, Celina MTM, Esper GK, Cláudio SP, José EL. Low sensitivity of NS1 protein tests evidenced during a dengue type 2 virus outbreak in Santos, Brazil, in 2010. Clinical and Vaccine Immunology 2012; 19(12): 1972-1976.
5. Maria GG, Thomas J, Roger G, Vo TTH, Shamala DS, Axel K, Susana V, Didyé R, Eric M, Juan CM, Angel B, Eva H, Efrén D, Prisca SAL, Sutee Y, Elci V, Herminia B, Iris V, Jeremy F, Cameron PS. Multi-country evaluation of the sensitivity and specificity of two commercially-available NS1 ELISA assays for dengue diagnosis. PLOS Neglected Tropical Diseases 2010; 4(8): e811, p1-10.
6. Guzman MG, Kouri G. Dengue diagnosis, advances and challenges. International Journal of Infection Diseases 2004; 8: 69-80.
7. TDR/WHO. In Dengue: Guidelines for diagnosis, treatment, prevention and control. TDR/WHO, Geneva, Switzerland, 2009.
8. PAHO. Dengue and Dengue Hemorrhagic Fever in the Americas: Guidelines for prevention and control. Pan American Health Organization Washington, DC, USA, 1994.
9. CDC. Laboratory Guidance and Diagnostic Testing. Centers for Disease Control and Prevention, CDC, USA, 2012. Retrieved from http://www.cdc.gov/dengue/clinicalLab/laboratory.html
10. Tan YY, Shamala DS, Seok MW, Firoz A, Anowar H, Bijon KS. Development of ASSURE® Dengue IgA Rapid Test for the detection of anti-dengue IgA from dengue infected patients. Journal of Global Infectious Diseases, 2011; 3(3): 233-240.

11. Firoz A, Huzzatul M, Meer TA, Ruhul A, Shamala DS, Seok MW, Tan YY, Bijon KS, Anowar MH. Evaluation of ASSURE® Dengue IgA Rapid Test using dengue-positive and dengue-negative samples. Diagnostic Microbiology and Infectious Disease, 2010; 68: 339-344.
12. Masaru N, Pan CY, Tsai WH, Chan CH, Sanae M, Tomohiko T, Lim CK, Harn MR, Ichiro K. Evaluation of immunoglobulin A – capture enzyme-linked immunosorbent assay for serodiagnosis of dengue virus infection. Dengue Bulletin, 2006; 30: 157-161.

LIMITED EXPRESSED WARRANTY DISCLAIMER

The manufacturer makes no express warranty other than that the test kit will function as an *in vitro* diagnostic assay within the specifications and limitations described in the product Instruction Manual when used in accordance with the instructions contained therein. The manufacturer disclaims any warranty expressed or implied, including such expressed or implied warranty with respect to merchantability, fitness for use or implied utility for any other purposes. The manufacturer is limited to either replacement of the product or refund of the purchase price of the product. The manufacturer shall not be liable to the purchaser or third parties for any damage, injury or economic loss howsoever caused by the product in the use or in the application thereof.

TECHNICAL PROBLEMS / COMPLAINTS

Should there be a technical problem / complaint, please do the following:

1. Note the kit lot number and the expiry date.
2. Retain the kits and the results that were obtained.
3. Contact the nearest MP Biomedicals office or your local distributor.



MP Biomedicals Asia Pacific Pte. Ltd.
2 Pioneer Place
Singapore 627885
Tel. + 65 6775 0008
Fax No.: + 65 6775 4536
Email: enquiry_ap@mpbio.com



MP Biomedicals Germany GmbH
Thüringer Straße 15
37269 Eschwege
Germany
Tel. (49) 5651 921 204
Fax No. : (49) 5651 921 181
Email : diagnostics@mpbio.com

Regional Office:
MP Biomedicals Germany GmbH
Thüringer Straße 15
37269 Eschwege
Germany
Tel No.: +49 5651 921 204
Fax No.: +49 5651 921 181
Email: diagnostics@mpbio.com