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1. Towbin H., Staehlin T. and Gordon J. Electrophoretic transfer of proteins from polyacrylamide gels to nitrocellulose sheets: procedure and some applications. Proc. Natl. Acad. Sci. U.S.A. 1976; 76: 4350-4354.
2. Poesz BJ., Ruscetti FW., Gazdar AF., Bonn PA., Minna JD., and Gallo RC. Detection and isolation of type C retrovirus particles from fresh and cultured lymphocytes of a patient with cutaneous T-cell lymphoma, Proc. Natl. Acad. Sci. U.S.A. 1980; 77(12): 7415-7419.
3. Kalyanaraman VS., Sarngadharan MG., Robert-Guroff M., Miyoshi I., Blayney D., Golde, and Gallo RC. A new subtype of human T-cell leukemia virus (HTLV-II) associated with a T-cell variant of hairy cell leukemia, Science 1982; 218: 571-573.
4. William AE., Fang CT., Slamon DJ. et al. Seroprevalence and epidemiological correlates of HTLV-I infection in U.S. blood donors, Science 1988; 240: 643-646.
5. Lee H., Swanson P., Shorty VS., Zack JA., Roseblatt JD. and Chen ISY. High rate of HTLV-II infection in seropositive IV drug abusers in New Orleans, Science 1989; 244: 471-475.
6. Lipka JJ., Bui K., Reyes GR., Moeckli R., Wiktorski SZ., Blattner WA., Murphy EL., Hanson CV., Shaw GM., Shinsky JJ., and Fong SKH. Determination of a unique immunodominant epitope of HTLV-I, Infect Dis 1990; 162: 353-357.
7. Wiktorski SZ., Alexandra SS., Shaw GM. et al. Distinguishing between HTLV-I and HTLV-II by Western Blot. Lancet 1990; 335: 1533.
8. Samuel KP., Lautenberger JA., Jorczyk CL., Josephs S., Wong Staal F. and Papas TS. Diagnostic potential for human malignancies of bacterially produced HTLV-I envelope protein. Science 1984; 226: 1094-1097.
9. Hadlock KG., Goh CJ., Bradshaw PA., Perkins S., Lo J., Habbaz RK., Kaplan J., and Fong SKH. Delineation of an immunodominant and highly HTLV specific epitope within the HTLV-I transmembrane glycoprotein, Blood 1995; 68(4): 1392-1399.

10. Varma M., Rudolph D., Knuchel M., Switzer W., Hadlock KG., Velligan M., Chan L., Fong SKH., Lai RB. Enhanced specificity of truncated transmembrane protein for serologic confirmation of HTLV-I and HTLV-II infection by Western Blot assay containing recombinant envelope glycoproteins, *J. Clin. Micro.* 1995; 33(12): 3239-3244.
11. Lillehoj EP., Alexander SS., Dubrule CJ., Wiktor S., Adams R., Thi A., Manns CC., and Blattner WA. Development and evaluation of a human T-Cell leukemia virus type I serologic confirmatory assay incorporating a recombinant envelope polypeptide, *J. Clin. Microbiol.* 1990; 28: 2653-2658.
12. Lal RB., Brodine SK., Coligan JE., and Roberts CR. Differential antibody responsiveness to p19 gag results in serological discrimination between human T-lymphotropic virus type I and type II, *J. Med. Virol.* 1991; 1: 232-236.
13. Hjelle B., Cyrus S., Swenson S., and Mills R. Serologic distinction between human T-lymphotropic virus (HTLV) type I and HTLV type II, *Transfusion* 1991; 31: 731-736.
14. Madeleine MM., Wiktor SZ., Goedert JJ., Manns A., Levine PH., Biggar RJ., Blattner WA. HTLV-I and HTLV-II worldwide distribution: reanalysis of 4,832 immunoblot results, *Int. J. Cancer* 1993; 54(2): 255-260.
15. De Castro-Costa CM, Araujo AQ, Barreto MM, Takayanagui OM, Sohler MP, da Silva EL, et al. Proposal for diagnostic criteria of tropical spastic paraparesis/HTLV-I-associated myelopathy (TSP/HAM), *AIDS Res Hum Retroviruses.* 2006;22(10):931-5.
16. Lal RB., Rudolph DL., Coligan JE., Brodine SK., and Roberts CR. Failure to detect evidence of human T-lymphotropic virus (HTLV) type I and type II in blood donors with isolated gag antibodies to HTLV-I/II, *Blood* 1992; 80: 544-550.
17. Khabbaz, RF., Heneine W., Grindon A., Hartley TM., Shulman G., and Kaplan J. Indeterminate HTLV serologic results in U.S. blood donors: Are they due to HTLV-I or HTLV-II? *Acquir. Immune Defic. Syndr.* 1992; 5: 400-404.
18. Lipka, JJ., Young KK., Kwok SY., Reyes GR., Snisky JJ, and Fong SK. Significance of human T-lymphotropic virus type I indeterminate serological findings among healthy individuals. *Vox Sang.* 1991; 61: 171-176.
19. Zrein M., Louwagie J., Boeykens H., Govers L., Hendickx G., Bosman F., Sablon E., Demarquilly C., Boniface M., and Saman E. Assessment of a new immunoassay for serological confirmation and discrimination of human T-cell lymphotropic virus infections. *Clin Diag. Lab. Immunol.* 1998;5: 45-49.
20. Witt DJ., Kuramoto K., Kemper M., and Holland P. Utility of prospective study of donors deferred as HTLV indeterminate. *Vox Sang.* 2000; 78:130-131.
21. Hayes C.G., Burans JP., and Oberst RB. Antibodies to Human T Lymphotropic virus Type I in a population from the Philippines: Evidence for cross-reactivity with *Plasmodium falciparum*. *The J. Infect. Dis.* 1990; 163: 257-262.
22. Gallo D., Diggs JL., and Hanson CV. Evaluation of two commercial Human T-Cell Lymphotropic Virus Western blot (Immunoblot) kits with problem specimens. *J. Clin. Microbiol.* 1994; 32: 2046-2049.
23. Garin B., Gosselin S., de The G., and Gessain A. HTLV-I/ II infection in a high viral endemic area of Zaire, Central Africa: Comparative evaluation of serology, PCR, and significance of indeterminate Western blot pattern. *J. Med. Virol.* 1994; 44: 104-109.
24. Fujiyama C., Fujiyoshi T., Matsumoto D., Yashiki S., Tamashiro H., and Sonoda S. Re-evaluation of anti-HTLV-I Western blot assay using HTLV-I and HTLV-II serum panels. *Clin. & Diag. Virol.* 1995; 4: 149-161.
25. Rouet F., Meertens L., Courouble G., Herrmann-Storck C., Pabingui R., Chancerel B., Abid A., Strobel M., Mauclere P., and Gessain A. Serological, epidemiological, and molecular differences between human T-cell lymphotropic virus type I - seropositive healthy carriers and persons with HTLV-I gag indeterminate Western blot patterns from the Caribbean. *J. Clin. Microbiol.* 2001; 39: 1247-1253.
26. Cesaire R., Bera O., Maier H., Lezin A., Martial J., Ouka M., Kerob-Bauchet B., Ould Amar AK., and Vernant JC. Seroindeterminate patterns and seroconversions to human T-lymphotropic virus type I positivity in blood donors from Martinique, French West Indies. *Transfusion* 1999; 39: 1145-1149.
27. Soldan SS., Graf MD., Waziri A., Flerlage AN., Robinson SM., Kawaninshi T, Leist TP., Lehky TJ., Levin MC., and Jacobson S. HTLV-I/II seroindeterminate western blot reactivity in a cohort of patients with neurological disease. *J. Infect. Dis.* 1999; 180: 685-694.
28. Taylor GP. The epidemiology of HTLV-I in Europe. *J Acquir Immune Defic Syndr Hum Retrovirol.* 1996;13 Suppl 1:S8-14.
29. World Health Organization. (2021). Human T-lympho-tropic virus type 1: technical report. World Health Organization. <https://iris.who.int/handle/10665/339773>. License: CC BY-NC-SA 3.0 IGO
30. Hino S, Yamaguchi K, Katamine S, et al. Mother-to-child transmission of human T-cell leukemia virus type-I. *Jpn J Cancer Res (Gann)* 1985;76:474-80.
31. Tajima K, Tominaga S, Suchi T, et al. Epidemiological analysis of the distribution of antibody to adult T-cell leukemia virus-associated antigen: possible horizontal transmission of adult T-cell leukemia virus. *Jpn J Cancer Res (Gann)* 1982; 73:893-901.
32. Okochi K, Sato H, Hinuma Y. A retrospective study on transmission of adult T cell leukemia virus by blood transfusion: seroconversion in recipients. *Vox Sang* 1984; 46:245-53.
33. Smadja D, Bellance R, Cabre P, Arfi S, Vernant JC. Clinical characteristics of HTLV-1 associated dermato-polymyositis. Seven cases from Martinique. *Acta Neurol Scand.* 1995;92(3):206-12.
34. Einsiedel L, Fernandes L, Spelman T, Steinfort D, Gotuzzo E. Bronchiectasis is associated with human T-lymphotropic virus 1 infection in an Indigenous Australian population. *Clin Infect Dis.* 2012;54(1):43-50.

35. Merle H, Cabre P, Olindo S, Merle S, Smadja D. Ocular lesions in 200 patients infected by the human T-cell lymphotropic virus type 1 in martinique (French West Indies). *Am J Ophthalmol.* 2002;134(2):190-5.
36. La Grenade L. HTLV-I, infective dermatitis, and tropical spastic paraparesis. *Molecular Neurobiology.* 1994; 8:147-53.

製造業者名稱：MP Biomedicals Asia Pacific Pte Ltd.
 製造業者地址：2 Pioneer Place, Singapore 627885
 聯絡號碼. : + 65 6775 0008
 傳真號碼. : + 65 6774 6146
 郵件地址：enquiry_ap@mpbio.com

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