



ALTERNATIVE PROCEDURE - OVERNIGHT ASSAY

Procedure:

1. Add 2 ml of DILUTED WASH BUFFER to each well.

2 ml
2. Using forceps, carefully remove required number of STRIPS from the tube and place numbered side up into each well. Include strips for Strong Reactive and Non-Reactive controls.

2 minutes
3. Incubate the strips for 1 to 2 minutes at room temperature (25 ± 3°C) on a rocking platform (speed of 12 to 16 cycles per minute). Remove buffer by aspiration. (Note: Do not allow the strips to dry. Failure may result in watery marks on developed strips for some specimens.)

2 minutes
4. Add 2 ml of BLOTTING BUFFER to each well.

2 ml
5. Add 20 µl each of patients' sera or controls to appropriate wells.

20 µl
6. Cover the tray with the cover provided and incubate overnight (16 - 20 hours) at room temperature (25 ± 3°C) on the rocking platform.

overnight
7. Carefully uncover the tray to avoid splashing or mixing of samples. Tilt the tray to aspirate the mixture from the wells. Change aspirator tips between samples to avoid cross-contamination.

3 x 2 ml
8. Wash each strip 3 times with 2ml of DILUTED WASH BUFFER allowing 5 minutes soak on the rocking platform between each wash.

3 x 2 ml
9. Add 2 ml of WORKING CONJUGATE SOLUTION to each well.

2 ml
10. Cover tray and incubate for 30 minutes at room temperature (25 ± 3°C) on the rocking platform.

30 minutes
11. Aspirate CONJUGATE from the wells. Wash as in step 8.

3 x 2 ml
12. Add 2 ml of SUBSTRATE SOLUTION to each well.

2 ml
13. Cover tray and incubate for 15 minutes on the rocking platform. (Note: The reaction can be stopped before 15 minutes if all the bands are visible.)

15 minutes
14. Aspirate the SUBSTRATE and rinse the strips at least three times with reagent grade water to stop the reaction (A dark background can result if washing is insufficient at this step).

3 x 2 ml
15. Using forceps, gently remove strips onto paper towels. Cover with paper towels and dry. Alternatively, allow strips to dry in the wells of the tray.

16. Mount strips on worksheet (non-absorbent white paper). Do not apply adhesive tape over the developed bands. Observe the bands (See Interpretation of Results) and grade the results. For storage, keep the strips in the dark.

SUMMARY OF ASSAY PROTOCOLS				
Reagents	Qty	Room Temp Rapid Assay	Room Temp Overnight Assay	
Nitrocellulose strip	1	-	-	
Wash Buffer	2 ml	1 - 2 mins	1 - 2 mins	
Blotting Buffer	2 ml	-	-	
Specimen	20 µl	60 mins	Overnight (16 - 20 hours)	
Wash Buffer	3 x 2 ml	3 x 5 mins	3 x 5 mins	
Conjugate	2 ml	60 mins	30 mins	
Wash Buffer	3 x 2 ml	3 x 5 mins	3 x 5 mins	
Substrate (Ready to use)	2 ml	15 mins (or less)	15 mins (or less)	
Distilled Water	3 x 2 ml	-	-	

AMOUNT OF REAGENTS REQUIRED FOR VARIOUS NUMBER OF STRIPS								
Reagents	NUMBER OF STRIPS TO BE USED							
	3	6	9	15	20	27	36	
1X Wash Buffer (ml)	60	100	140	240	300	400	520	
1X Blotting Buffer (ml)	20	40	60	80	100	120	160	
Conjugate (µl)	11	17	23	35	45	59	77	
Substrate (ml)	11	17	23	35	45	59	77	
Blotting Powder (g)	1	2	3	4	5	6	8	

QUALITY CONTROL

We recommend that the Non-Reactive & Strong Reactive controls be run with every batch of assay. In order for the results obtained from any assay to be considered valid, the following conditions must be met:

1. **NON-REACTIVE CONTROL**

No viral specific bands should be observed on the Non-Reactive control strip.
2. **STRONG REACTIVE CONTROL**

All relevant molecular weight bands must be evident. Figure 1 provides a guide to the relative positioning of bands visualized with the MP Diagnostics HIV-2 BLOT 1.2 kit and permits identification of bands observed for the STRONG REACTIVE CONTROL. The bands are p26, gp36, p53, p56, p68, gp80 and gp125.

INTERPRETATION OF RESULTS

The major HIV-2 gene products that have been identified are:

GAG	POL	ENV
CORE ANTIGENS	ENDONUCLEASE POLYMERASE ANTIGENS	ENVELOPE ANTIGENS
p16, p26, p56	p34, p53, p68	gp36, gp80, gp125

The criteria for interpretation of a positive HIV-2 specimen may differ according to various national regulatory bodies. It is recommended that national policies be followed in the interpretation of results. However, the recently proposed WHO criteria for interpretation has been widely used. The WHO criteria can be used in the interpretation of results of test specimens performed with the MP Diagnostics HIV-2 BLOT 1.2.

Each strip is compared to the strips used with the NONREACTIVE AND STRONG REACTIVE CONTROLS for that run.

The result of blotting is then interpreted as NEGATIVE, INDETERMINATE or POSITIVE based on the pattern which is present:

PATTERN	INTERPRETATION
1) No viral specific bands present	NEGATIVE
2) 2 of the 3 envelope bands (gp36, gp80, gp125)	POSITIVE
3) Any bands present but pattern does not meet criteria for POSITIVE	INDETERMINATE

Individuals seropositive for either HIV-1 or HIV-2 may display cross-reactivity with either viral proteins. The GAG and POL genes of HIV-1 and HIV-2 are highly conserved, resulting in a high degree of serological cross-reactivity with the GAG and POL encoded proteins. There is a greater divergence for the amino acids of the envelope glycoproteins and cross-reactivity against those glycoproteins are considered to occur less often.

However reports (4, 5) have shown that cross-reactivity against the envelope glycoproteins occur more often than previously believed. Thus sera showing extensive cross-reactivity indicating dual infection should be interpreted carefully. Further confirmation by characterization of gene sequences via polymerase chain reactions (PCR) procedure is recommended in such highly cross-reactive specimens to indicate infection with either HIV-1, HIV-2 or dually infected.

LIMITATIONS OF THE PROCEDURE

Optimal assay performance requires the strict adherence to the assay procedure described. Deviation from the procedure may lead to aberrant results.

Individuals with POSITIVE blots for antibodies to HIV-2 should be referred for medical evaluation. A NEGATIVE blot does not exclude the possibility of infection with HIV-2. INDETERMINATE blots should not be used as the basis for diagnosis of HIV-1 infection. It is recommended that REPEAT TESTING using the original specimen or FOLLOW-UP TESTING be performed for such sera.

LIMITED EXPRESSED WARRANTY DISCLAIMER

The manufacturer makes no express warranty other than that the test kit will function as an *in-vitro* diagnostic assay within the specifications and limitations described in the product Instructions For Use when used in accordance with the instructions contained therein. The manufacturer disclaims any warranty expressed or implied, including such expressed or implied warranty with respect to merchantability, fitness for use or implied utility for any other purposes. The manufacturer is limited to either replacement of the product or refund of the purchase price of the product. The manufacturer shall not be liable to the purchaser or third parties for any damage, injury or economic loss howsoever caused by the product in the use or in the application thereof.

TECHNICAL PROBLEMS / COMPLAINTS

Should there be a technical problem/complaint, please do the following:

1. Note the kit lot number and the expiry date.

2. Retain the kits and the results that were obtained.

3. Contact the nearest MP Biomedicals office or your local distributor.

REFERENCES

1. WHO Collaborating Group on HIV-2, WHO Weekly Epidem Rec 1990; 10:74-75.

2. Clavel F., Guetard D., Brun-Vezinet F., et al. Isolation of a new human retrovirus from West African patients with AIDS. Science 1986; 223:343-346.

3. Clavel F., HIV-2, the West African AIDS virus. AIDS 1987; 1:135-140.

4. R.S. Tedder, Hughes A., Corrah T., et al. Envelope cross-reactivity in Western Blot for HIV-1 and HIV-2 may not indicate dual infection. Lancet 1988; 11:927-930.

5. Bottiger B., Karlsson A., Andreasson F., et al., Envelope cross-reactivity between Human Immunodeficiency Virus Type 1 and Type 2 detected by Different Serological Methods: Cross Reactivity against the main neutralizing site. J. Virol. 1990; 65(7):3492-3499.



**MP Biomedicals Asia Pacific Pte Ltd**  
2 Pioneer Place  
Singapore 627885

Tel No. : + 65 6775 0008  
Fax No. : + 65 6774 6146  
Email : enquiry\_ap@mpbio.com

Regional Office:

**MP Biomedicals Germany GmbH**  
Thüringer Straße 15  
37269 Eschwege  
Germany  
Tel. No.: +49 5651 921 204  
Fax No.: +49 5651 921 181  
Email: diagnostics@mpbio.com

FIGURE 1

