

MagBeads FastDNA[®] Kit

Ready-to-Use for MPure-32™ aNAP System



Cat. No.: 117033600 (96 PREPS)

Quick-Start Manual

Revision 3.0 Dec 2023

Notes before starting:

- Add 2.5 mL Protease Dissolve Buffer to the Proteinase K, and store at -20~8°C after dissolve.
- Add 1.25 mL Protease Dissolve Buffer to the RNase A, and store at -20~8°C after dissolve.

Automation Extraction

Sample Preparation

A. Solid tissue (1-20mg)

1. Cut ~20 mg tissue into small pieces and transfer into a new 1.5 mL microcentrifuge tube. Add 200 μ L Buffer ATL and 20 μ L Proteinase K and incubate with shaker at 55°C for 30 to 180 mins.
2. (Optional) Add 10 μ L RNase A into the lysate and incubate at room temperature for 10 mins.
3. Add 200 μ L Buffer AL into the samples, mix thoroughly by vortexing to yield homogenous solution.
4. Incubate at 70°C for 10 mins.

B. Anticoagulated blood or Plasma (200 μ L)

1. Transfer 20 μ L Proteinase K to a new 1.5 mL microcentrifuge tube.
2. Add 200 μ L whole blood, plasma, or other body fluids into the tube. Shake to mix for 5 seconds.
3. Add 200 μ L Buffer AL into the samples. Invert the tube for 3 - 5 times, and vortex at maximum speed for 10s. Incubate at 70°C for 10 mins.

C. Saliva sample

1. Add 20 μ L Proteinase K and 10 μ L RNase A into a 1.5 mL microcentrifuge tube.
2. Transfer 450 μ L saliva to the tube and shake to mix for 5 s.
3. Incubate at 55°C for 30 mins.

D. Culture cells

1. Collect cells (<2 x10⁶) by centrifuging at 2,000 x g for 5 mins. Remove the supernatant.
2. Add 200 μ L Buffer PBS, 20 μ L Proteinase K, and 10 μ L RNase A into the sample. Resuspend the cells by vortexing.
3. Add 200 μ L Buffer AL and vortex for 10 seconds. Incubate the mixture at 70°C for 10 mins.

E. Semen sample

1. Transfer 100 μ L semen to a 1.5 mL microcentrifuge tube.
2. Add 100 μ L Buffer ATL, 10 μ L DTT Solution (1M), and 20 μ L Proteinase K into the samples. Shake at 55°C for 30 mins.
3. Add 200 μ L Buffer AL into the sample, then vortex to mix and incubate at 70°C for 10 mins.

F. Swab DNA extraction

1. Transfer the swabs to a 2.0 mL microcentrifuge tube.
2. Add 500 μ L ATL and 20 μ L Proteinase K to the swab. Shake at 55°C for 15 to 30 mins.
3. Transfer the supernatant into a new 2.0 mL microcentrifuge tube.

G. Blood stains/Seminal stains

1. Transfer 3 slices (3mm) to the 2.0 mL microcentrifuge tube.
2. Add 250 μ L Buffer ATL and 20 μ L Proteinase K to the sample. Shake at high speed for 30 to 60 mins at 55°C.
3. Add 250 μ L Buffer AL into the samples, Shake at high speed for 10 mins at 70°C.
4. Centrifuge at 13,000 x g for 1 min.
5. Transfer 400 μ L of the supernatant to a new 2.0 mL microcentrifuge tube.

H. FFPE Samples

1. Using a scalpel, trim excess paraffin off the sample block. Transfer 1 to 3 sections (5 - 10 μ m) into a 1.5 mL microcentrifuge tube.
2. Remove Paraffin by xylene or Buffer DPS (not provided).
3. Add 200 μ L Buffer ATL and 20 μ L Proteinase K into the sample, mix well and incubate at 56°C for 60 mins, and 90°C for 60 mins.
4. Cool the sample to room temperature and add 200 μ L Buffer AL. Mix well before proceeding to extraction step.

MagBeads FastDNA[®] Kit

Ready-to-Use for MPure-32™ aNAP System



Cat. No.: 117033600 (96 PREPS)

Quick-Start Manual

Revision 3.0 Dec 2023

Extraction

1. Transfer 400 µL of supernatant carefully to well #1 or #7 of the pre-filled reagent.
2. Place the reagent plate on MPure-32™ aNAP System and run the assay with the program named "FastDNA" which has the following setting:

Step	Well	Process	Time (s)			Mixing Speed	Temp (°C)
			Mix	Wait	Attract		
1	#1/#7	Bind	240	0	90	Medium	RT
2	#2/#8	Wash 1	120	0	60	Medium	RT
3	#3/#9	Wash 1	90	0	60	Medium	RT
4	#4/#10	Wash 2	90	0	60	Medium	RT
5	#5/#11	Wash 2	90	0	60	Medium	RT
6	#5/#11	Dry	0	300	0	-	RT
7	#6/#12	Elute	480	0	120	Medium	55
8	#1/#7	Magbeads Release	60	0	0	Medium	RT

3. Transfer the eluted DNA into a clean 1.5 mL microcentrifuge tube. DNA is now ready for PCR and other downstream applications. Store the purified nucleic acid at -20°C for an extended storage.

Note: If there are still Magnetic Beads remaining in eluted DNA, please centrifuge at 14,000 x g for 3-5 mins and transfer the supernatant into a clean 1.5 mL microcentrifuge tube.

Ordering Information

Equipment

Automated extraction system from low to high throughput

Catalog No.	Product Name	Throughput
07EMC043	MPure-32™ aNAP System	Up to 32 samples
07EMC044	MPure-96™ aNAP System	Up to 96 samples



Instruments for lysing and homogenizing environmental and biological samples

Catalog No.	Product Name
116004500	FastPrep™ Classic
116005500	FastPrep-24™ 5G
116010500	FastPrep-96™
116012500	SuperFastPrep-2™



MP BIOMEDICALS

APAC: +65 6775 0008 | custserv.ap@mpbio.com
 EUROPE: 00800 777 9999 | custserv.eur@mpbio.com
 AMERICAS: 800 854 0530 | custserv.na@mpbio.com
 Learn more at www.mpbio.com



Scan for detailed instruction manual

