

MP Biomedicals

Spin-Column Extraction Kit for Soil gDNA

SPINeasy DNA Kit for Soil



Soil is one of the most commonly studied environmental sample. It has great variety and rich in biological content. However, the process to extract DNA from soil is not an easy task. Soil particles hinder effective lysis of cell bodies, while humic substances in soil bind irreversibly to DNA. All of which will reduce the yield and purity of extracted DNA.

With decades of experiences in environmental handling, we dedicate to produce the most effective tools for fellow researchers. Proudly presenting the SPINeasy DNA Kit for Soil, a kit specialized to extract gDNA of high yield and purity from various soil types.



Product Highlights

01

Lysing Matrix E tube with three kinds of glass beads can effectively lyse cell bodies to release DNA, ensuring non-biased study of the microbiome.

02

Proprietary lysis buffers to protect released DNA and remove unwanted RNA.

03

Specially formulated buffers to get rid of humic substances and other inhibitors, greatly improve the purity of extracted gDNA.


04

Column S1 with high binding capacity and selectivity for gDNA. It is designed with large volume to reduce frequency of liquid transfer.



- 1 Apply for trial kit at our official website.
- 2 Feedback after trial and receive a complimentary gift.

Ordering Information

Product Name	Format	Cat. No.	Method	Characteristic
 SPINeasy DNA Kit for Soil	50 preps	116530050	Spin-Column	Suitable for various soil types, easy to operate
SPINeasy DNA Kit for Soil (Sample)	5 preps	116530000		
FastDNA Spin Kit for Soil	50 preps	116560200	Glassmilk	Suitable for soil type with extremely low biomass content
MagBeads FastDNA Kit for Soil	50 preps	NA	Magnetic Beads	Upcoming product in Sept 2020. Suitable for various soil types, automation possible



MP Biomedicals

Spin-Column Extraction Kit for Soil gDNA

1 Quality and Quantity of Extracted gDNA

Table 1: Quality and quantity of gDNA extracted from various soil samples using SPINeasy DNA Kit for Soil.

Sample Type	Yield(ng/mg sample)	260/280	260/230
Organic soil	47.78±0.21	1.79±0.01	1.10±0.02
Paddy soil	16.03±0.06	1.94±0.01	1.70±0.04
Flowerbed soil	14.88±0.65	1.84±0.01	1.20±0.04
Saline soil	7.55±0.65	1.89±0.01	1.40±0.00
Desert soil	1.64±0.19	1.85±0.06	1.02±0.00

SPINeasy DNA Kit for Soil is compatible with various soil types. Extracted gDNA is of high yield and purity.

2 Electrophoretic Analysis




Figure 1: gDNA extracted from different soil samples using SPINeasy DNA Kit for Soil
 M: 1kb plus DNA ladder
 Lane 1: Organic soil
 Lane 2: Paddy soil
 Lane 3: Flowerbed soil
 Lane 4: Saline soil
 Lane 5: Desert soil

Extracted gDNA is highly intact.

3 16S rDNA PCR amplification






Figure 2: PCR amplification of 16S rDNA gene from different soil samples using SPINeasy DNA Kit for Soil


M: 1kb plus DNA ladder
 Lane 1: Organic soil
 Lane 2: Paddy soil
 Lane 3: Flowerbed soil
 Lane 4: Saline soil
 Lane 5: Desert soil
 Lane 6: Negative control

Extracted gDNA is free from inhibitors and ready for downstream analysis.

4 Performance Comparison



Remark 1: Pure DNA has 260/280 at 1.80 (Indicated with red line).



Remark 2: Pure DNA has 260/230 > 1.00.

SPINeasy DNA Kit for Soil has consistently achieved better yield and purity as compared to competitor kits.

MP Bio Americas

custserv@mpbio.com
 1.800.854.0530

MP Bio Europe

custserv.eur@mpbio.com
 00800.7777.9999

MP Bio APAC

asia.custserv@mpbio.com
 65.6775.0008

MP
 www.mpbio.com