


AutoScan Pro 4.0



Experience Automated and Accurate Western and Line Blot Test Analysis & Interpretation with our Latest Software

Accurate result interpretation
in 4 simple steps

Objective and Standardized
interpretation of Western Blot and
Immunoblot results

Easy integration with your
current Laboratory Information
Management System

MP Bio Americas

📧 custserv@mpbio.com
📞 1.800.854.0530

MP Bio Europe

📧 custserv.eur@mpbio.com
📞 00800.7777.9999

MP Bio APAC

📧 asia.custserv@mpbio.com
📞 65.6775.0008




AutoScan Pro 4.0

AutoScan Pro 4.0 is a software which enables an automated analysis and interpretation of Western Blot and Line Immunoblot test results. AutoScan Pro 4.0 reads the digital image of the strips generated by a validated scanner. Using the pre-programmed interpretation criteria recommended by MP Biomedicals, AutoScan Pro 4.0 evaluates the test strip results and provides an accurate interpretation to the user.

- Allows full use of the lot-specific protein finder which is embedded in a barcode. Every barcode is unique and specific to a kit lot, hence enabling you to get the most accurate analysis and result interpretation for every test.
- Enables an objective reading and interpretation of Western Blot and Line Immunoblot results.
- Highly sensitive and precise. Bands of different intensities are smartly detected and individually compared with the control strip.
- Includes comprehensive patient and result search functions, which helps in retrieving and analyzing the results from a specific patient cohort.
- Comes with free extended database management functions which allow you to analyze, interpret and generate the patient reports with a single software.

Accurate result interpretation in 4 easy steps:

- 1 Paste the strips on the analysis template
- 2 Scan the lot-specific barcode
- 3 Scan the analysis template and run analysis with AutoScan Pro 4.0
- 4 View and export your results conveniently



Patient Report & Optional Data Interface to LIS



Patient	Patient-ID	Sample-ID	Test	Image	Result
Karsten, Christian	9494	00182	MP Biomedicals HIV BLOT 2.2 9A 00182		HIV-1 Positive with HIV-2 indicated
Zahn, Eduardo	NCE	00183	MP Biomedicals HIV BLOT 2.2 9A 00182		HIV-1 Positive with HIV-2 indicated
Junker, Sabine	7541	00184	MP Biomedicals HIV BLOT 2.2 9A 00182		HIV-1 Positive with HIV-2 indicated
Euda, Alo	422631	00185	MP Biomedicals HIV BLOT 2.2 9A 00182		HIV-1 Positive with HIV-2 indicated

Description	Cat. No.
AutoScan Pro 4.0	07EMC021
AutoScan Pro 4.0 (Upgrade from 3.4)	07EMC022

LEARN MORE
www.mpbio.com

A special upgrade offer is available to existing users of AutoScan Pro 3.4. Please contact a sales representative or distributor in your region to find out more.

MP Bio Americas
custserv@mpbio.com
1.800.854.0530

MP Bio Europe
custserv.eur@mpbio.com
00800.7777.9999

MP Bio APAC
asia.custserv@mpbio.com
65.6775.0008

