

Simplify Your Growth Media Essentials

- ✓ **EASY-TO-USE**
Simply drop the capsules in water and autoclave
- ✓ **EXCEPTIONAL HIGH CELL DENSITIES**
and plasmid yields
- ✓ **WIDE SELECTION OF MEDIA FORMULATIONS**
both agar and broth formulations available
- ✓ **MOLECULAR BIOLOGY CERTIFIED™**
for the highest quality



SELECT YOUR MEDIA TYPE:

Bacterial Growth Media	Bacterial media for higher density growth and maintenance of recombinant <i>E. coli</i> strains. ▶ CircleGrow, LB, LB-agar, L-Broth, L-Broth agar, Terrific Broth I and II, SOB
Lambda and Filamentous Bacteriophage Media	Optimal media for the cultivation of recombinant <i>E. coli</i> and propagation of λ bacteriophage. ▶ SuperBroth, NZCYM Broth, YT Broth, YT-agar Broth, 2XYT Broth
Yeast Growth Media	Commonly used yeast media for maintenance and propagation of <i>P. pastoris</i> and <i>S. cerevisiae</i> yeast strains. ▶ YPD Broth, YPD-agar Broth, YEP Broth, SD Medium
Other Media	A selection of media to meet your needs for bacteria, yeast, and phage media. ▶ CSM, Tryptic soy, Columbia agar, MacConkey agar



WANT TO LEARN MORE?

Scan the code to download a digital copy of this flyer for a complete growth media product listing



Learn more at www.mpbio.com

Product Name	Format	Size	Cat. No.
Bacterial Growth Media			
CIRCLEGROW BROTH 40 g/L; Content: <i>proprietary formulation</i>	Capsules	250 or 1000 capsules	113000101
	Large Capsules	227 g / 454 g / 1 kg	113000116
	Pouch	0.5 L or 1 L pouches (10 ea)	113000165
	Media	500 mL	113000104
LB Medium 25 g/L; Content per liter: 10 g tryptone, 5 g yeast extract, 10 g NaCl	Capsules	227 g / 454 g / 1 kg / 5 kg	113002011
	Large Capsules	227 g / 454 g / 1 kg / 5 kg / 15 kg	113002016
	Pouch	0.5 L or 1 L pouches (10 ea / 100 ea)	113002065
	Media	500 mL	113002014
LB-agar Medium 40 g/L; Content per liter: LB Medium, 15 g Agar	Capsules	227 g / 454 g / 1 kg / 5 kg	113002201
	Large Capsules	227 g / 454 g / 1 kg / 5 kg	113002216
	Pouch	0.5 L or 1 L pouches (10 ea)	113002265
LB Medium, Lennox 20 g/L; Content per liter: 10 g tryptone, 5 g yeast extract, 5 g NaCl	Capsules	454 g / 1 kg / 5 kg	113002121
	Large Capsules	1 kg / 5 kg	113002136
LB Agar, Lennox 35 g/L; Content per liter: LB Medium Lennox, 15 g Agar	Capsules	454 g / 1 kg	113002321
	Large Capsules	454 g	113002326
L-Broth 5.5 g/L; Content per liter: 10 g tryptone, 5 g yeast extract, 0.5 g NaCl	Capsules	1 kg / 5 kg	113001021
L-Broth Agar 30.5 g/L; Content per liter: L-Broth, 15 g Agar	Capsules	454 g / 1 kg / 5 kg	113001221
	Large Capsules	1 kg	113001236
Terrific Broth Part 1: 36 g/L; Part 2: potassium phosphate salts and glycerol in sterile solution; Content per liter: 24 g yeast extract, 12 g tryptone	Capsules	1 kg	113045031
Terrific Broth II 50 g/L; Content: <i>proprietary formulation</i>	Capsules	454 g / 1 kg / 5 kg	113046021
	Pouch	1 L pouches (10 ea)	113046075
SOB Medium 30.5 g/L; Content per liter: 20 g tryptone, 5 g yeast extract, 0.5 g NaCl, 5 g MgSO ₄ •7H ₂ O	Capsules	454 g	113030021

Product Name	Format	Size	Cat. No.
Lambda and Filamentous Bacteriophage Media			
SuperBroth 60 g/L; Content per liter: 35 g tryptone, 20 g yeast extract, 5 g NaCl	Capsules	1 kg	113010031
	Large Capsules	1 kg	113010036
NZCYM Broth 23 g/L; Content per liter: 10 g NZ-Amine, 1 g casamino acids, 5 g yeast extract, 2 g MgSO ₄ •7H ₂ O, 0.5 g NaCl	Capsules	1 kg	113004031
YT-Broth 18 g/L; Content per liter: 8 g tryptone, 5 g yeast extract, 5 g NaCl	Large Capsules	5 kg	113011046
YT-Broth Agar 33 g/L; Content per liter: YT Broth, 15 g Agar	Capsules	1 kg / 5 kg	113011242
2xYT-Broth 31 g/L; Content per liter: 16 g tryptone, 10 g yeast extract, 5 g NaCl	Capsules	454 g / 1 kg / 5 kg	113012021
	Large Capsules	1 kg / 5 kg	113012037
Yeast Growth Media			
YPD (YEPD) Broth 50 g/L; Content per liter: 20 g peptone, 10 g yeast extract, 20 g dextrose	Large Capsules	1 kg	114001036
	Pouch	0.5 L or 1 L pouches (10 ea)	114001065
	Media	1 L	114001034
YPD (YEPD) Broth agar Content per liter: 20 g peptone Y; 10 g yeast extract-Y; 20 g dextrose; 17 g Agar-Y	Pouch	0.5 L or 1 L pouches (10 ea)	114001265
YEP Broth 30 g/L; Content per liter: 20 g peptone, 10 g yeast extract	Large Capsules	5 kg	114004046
SD medium-LEU	Pouch	0.5 L or 1 L pouches (10 ea)	114811065
SD medium-LEU-TRP	Pouch	0.5 L or 1 L pouches (10 ea)	114823065
SD medium-ADE-HIS-LEU-TRP	Pouch	0.5 L or 1 L pouches (10 ea)	114842065
SD medium-HIS-LEU-TRP	Pouch	0.5 L or 1 L pouches (10 ea)	114830065
SD medium-URA	Pouch	0.5 L pouches (10 ea)	114813065
SDA medium-LEU	Pouch	0.5 L pouches (10 ea)	114811165
SD Agar medium-ADE-HIS-LEU-TRP	Pouch	0.5 L pouches (10 ea)	114842165



DEPENDABILITY DELIVERED



ISO 13485: 2016



ISO 9001:2015



WHO GMP



Media



Buffer



Antibiotics
Solution



Custom Solution/
Reagents



Rapid
Tests

MP Biomedicals India Pvt Ltd,

Plot No C-57/3, TTC Industrial Area, Pawane Village, Turbhe, Navi Mumbai 400705.

☎ **+91.2227636921/22/25** | ✉ **custserv.in@mpbio.com** | **www.mpbio.com**