

Selection Guide MagBeads FastDNA & FastRNA Extraction and Purification Kits

Kit Name Sample Type		DNA Extraction							RNA Extraction					All Purpose		Other	
		FastDNA Kit	Blood	Soil	Feces	MaxPure Feces	FFPE	Circulating	Microbiome	FastRNA Kit	Blood	Feces	FFPE	Virus	Pathogen	Virus	Plasmid Miniprep
Animal Tissues, Clinical Samples & Cultured Cells	Fresh/Frozen Tissue	•								•							
	FFPE Tissue	•											•				
	Swabs	•								•				•			
	Cultured Cells	•							•							•	
	Cell Culture Media													•			
	Lymphocytes	•									•						
Blood & Other Fluids	Whole Blood	•	•						•	•				•	•		
	Serum		•						•	•				•	•		
	Plasma		•						•	•				•	•		
	Buffy Coat		•						•	•				•	•		
	Bone Marrow		•						•	•				•			
	Semen	•	•						•	•					•		
	Saliva	•							•	•				•	•		
	Soaking Solution/PBS		•						•	•				•	•		
Environmental & Stool Samples	Sediment			•					•								
	Sand			•					•								
	Soil			•					•								
	Feces			•	•	•			•			•					
	Intestinal Content	•		•	•	•											
	DSS-treated/inhibitor-rich feces					•				•		•					
Plants & Microbes	Plant Tissues	•							•	•							
	Fungi			•	•				•	•				•	•		
	Bacteria			•	•				•	•				•	•		•
	Protists			•	•									•	•		
	Viruses			•	•									•	•		
	Catalog No.*	116575192	116574096	116561050	116570400	116598050	116576096	116577192	116583000	116572096	116586050	116588050	116573192	116578050	116589050	116571192	116587050

*Catalog numbers listed refer to the manual bottled version; prefilled (ready-to-use) formats for automated nucleic acid systems are also available.