

Rapid SARS-CoV-2 Antigen Test Card



The Rapid SARS-CoV-2 Antigen Test Card is a rapid, chromatographic immunoassay for the qualitative detection of SARS-CoV-2 antigens present on human nasal, nasopharyngeal or oropharyngeal swab specimens.

- ▶ **EASY TO USE:** Easy sampling using nasal, nasopharyngeal or oropharyngeal swab specimens.
- ▶ **RAPID AND CONVENIENT:** Results in 15 minutes. Requires no instruments.
- ▶ **EXCELLENT PERFORMANCE CHARACTERISTICS:** Specificity 99.16%, Sensitivity 96.17%, Accuracy 98.25% (cumulated)



| | |
|---------------------|--|
| Package | 20 tests includes all required consumable materials including sterile swabs |
| Test Type | Qualitative |
| Instrument | No instrument required (Point of Care Test) |
| Specimens | Nasal, nasopharyngeal and oropharyngeal |
| Duration | Results after 15 minutes |
| Specificity | 99.14 % (nasopharyngeal) 99.20 % (oropharyngeal) 99.20 % (nasal) |
| Sensitivity | 96.08 % (nasopharyngeal) 96.23 % (oropharyngeal) 96.19 % (nasal) |
| Storage temperature | 2 - 30°C |
| Usage | For professional <i>in-vitro</i> diagnostic use only. |

This product is CE-IVD marked for professional use and complies with the European In-Vitro Diagnostic Directive (98/79/EC).

SHOP NOW

SKU: [07AG6020B](#)



MP Biomedicals Germany GmbH
 ✉ diagnostics@mpbio.com
 ☎ 0049.5651.921186

Learn more at www.mpbio.com

