



MP Biomedicals, LLC

29525 Fountain Parkway
Solon, Ohio 44139

Telephone: 440/337-1200

Toll Free: 800/854-0530

Fax: 440/337-1180

mailto: biotech@mpbio.com

web: <http://www.mpbio.com>

TECHNICAL INFORMATION

Catalog Number: 100921, 100924, 102566

Papain

Molecular Weight: 23,000

CAS # : 9001-73-4

Physical Description: White to off white powder

Description: Papain is a sulfhydryl protease from *Carica papaya* latex. (A second protease, chymopapain, and a lysozyme have also been isolated from this same source.) Since native crystalline papain is quite unreactive until acted upon by mild reducing agents such as cysteine or cyanide, it may exist as a zymogen (Brocklehurst and Kierstan 1973). For a general review, see Liener (1974). Barrett and Buttle (1985), Polgár (1984) and Brocklehurst and Salih (1983) report on the classification of papaya latex proteinases. Papain has wide specificity. In her review, Arnon (1970) has indicated that it will degrade most protein substrates more extensively than the pancreatic proteases. It is also an esterase. Papain has been reviewed by Smith and Kimmel (1960). It has been reported by Sluyterman and Wijdenes (1972a) that the action of papain on leucine methyl ester produces an insoluble polyleucine peptide. The finding of Thomas (1956) that papain breaks down the intercellular matrix of cartilage (see also McCluskey and Thomas 1958), led to its further study as a chondromucoproteinase (Smith et al. 1962). Proteolytic enzymes are widely used in cell isolation. With some tissues papain has proved less damaging and more effective than other proteases. Lam (1972) found that of the enzymes used for dissociating turtle retina, papain produced the least trauma. Intact single photoreceptor cells have also been isolated from adult salamander retina with papain (Bader et al. 1978, Townes-Anderson et al. 1985). Huettner and Baughman (1986) described a method using papain to obtain high yields of viable, morphologically intact cortical neurons from postnatal rats. Finkbeiner and Stevens (1988) applied the Huettner and Baughman method to the dissociation of postnatal rat hippocampus. Papain is used with fetal as well as postnatal brain regions to provide maximal dissociation and viability of neurons.

Composition: Papain is a single peptide chain of 211 residues folded into two parts that form a cleft (Dreuth et al. 1968). A three-dimensional structure has been indicated by Wolthers et al. (1970). The molecule has one free SH group which is functional (Smith et al. 1975; Shipton et al. 1975). According to Alecio et al. (1974) there are seven subsites each capable of accommodating a single amino acid residue of a peptide substrate. See also Glick and Brubacker (1974). Other reports on molecular information and its relation to activity are as follows: Fink and Gwyn (1974), Lewis and Shafer (1974), Akalski et al. (1973), Allen and Lowe (1973), Brocklehurst and Little (1973), Mole and Horton (1973a, b, and c), Banks and Shafer (1972), Brocklehurst et al. (1972), Campbell and Kaiser (1971, 1972), Sluyterman and Wijdenes (1972b), Hinkle and Kirsch (1971a and b) Jori et al. (1971), Lowe and Yuthavong (1971a and b) and Steiner (1971).

Optimum pH: 6.0 - 7.0.

Extinction coefficient: 25.0 (Mitchel et al. 1970).

Isoelectric point: pH 9.6 (Sluyterman and DeGraff 1972).

Activators: Papain is activated by cysteine, sulfide, sulfite, etc. It is enhanced when heavy metal binding agents such as EDTA are also present. Kirschenbaum (1971) indicated that N-bromosuccinimide enhances the activity. Hall et al. (1972) report on the affects of acridine dyes.

Inhibitors: Substances which react with sulfhydryl groups including heavy metals, carbonyl reagents (Moriyama 1967). Westerik and Wolfenden (1972) have studied aldehydes as papain inhibitors and Sluyterman and Wijdenes (1973) report on benzoylamidoacetonitrile as an inhibitor. See Shapira and Arnon (1967a and b) on antibody inhibitors. Papain may be inactivated by H₂O₂ generated by [γ]-irradiation of H₂O- the active SH group being oxidized to sulfenic acid. (Lin et al 1975). See also Allison and Swain (1973). Specific inhibitors are AEBSF, antipain, cystatin, E-64, leupeptin, PMSF, TLCK and TPCK.

Stability: Papain as a crystalline suspension is stable at 5°C for 6-12 months. Stabilizing agents are EDTA, cysteine and dimercaptopropanol.

To enhance stability as well as solubility it may be advantageous to convert crystalline papain to its mercury derivative (Brubacher and Bender 1966).

Assay

Method: A titrimetric determination of the acid produced during the hydrolysis of benzoyl-L-arginine ethyl ester (BAEE). One unit will hydrolyze one micromole of benzoyl-L-arginine ethyl ester per minute at 25°C and pH 6.2 under the specified conditions.

Reagents

Enzyme diluent (Activation buffer): Prepare fresh daily by mixing the following:

0.01 M EDTA	10 ml
0.06 M b-Mercaptoethanol	0.1 ml
0.05 M Cysteine-HCl	10 ml
Reagent grade water	70 ml

Substrate solution: Prepare fresh daily by mixing the following:

0.058 M BAEE	15.0 ml
0.01 M EDTA	0.8 ml
0.05 M Cysteine-HCl	0.8 ml

Adjust pH to 6.2 and dilute to a final volume of 21 ml with reagent grade water.

Titrant: 0.01-0.02 N NaOH, standardized

Enzyme

Activate enzyme by dissolving in enzyme diluent to a concentration of 0.05-0.1 mg/ml. Under these conditions activation is complete within 30 minutes.

Determination of protein concentration

Papain: mg protein/ml = $A_{280} \times 0.4$

Procedure

The reaction can be measured with either an automatic titrator or a laboratory pH meter. The titration vessel should be maintained at 25°C.

Pipette the following into the titration vessel at 25°C:

Substrate solution	5.0 ml
3.0 M NaCl	5.0 ml
Reagent grade water	5.0 ml

At zero time add 0.1 ml of appropriately diluted enzyme and adjust the pH to 6.2. Record the amount of standardized NaOH added per minute to maintain the pH at 6.2 after a constant rate is achieved.

Calculation

ml base added/minute x normality x 1000

U/mg = _

mg enzyme in reaction buffer

Technical note: Mercuripapain must be activated before use. Mercury is removed from the enzyme in activation buffer. After 30 minutes in this solution, the enzyme is completely activated and the mercury has been chelated. The mercuripapain suspension contains no free mercury. The product has been extensively dialyzed prior to packaging.

Availability:

Catalog Number

Description

Sizes

100921	Papain, 2X crystallized, suspension in 0.05 M sodium acetate. Activity is approximately 10 - 20 units/mg protein	25 mg 50 mg 100 mg 500 mg 1 g 5 g
100924	Papain, 2X crystallized, A suspension in 0.05 M sodium acetate. Activity is approximately 15 - 40 units/mg protein	25 mg 100 mg 250 mg 500 mg 1 g 5 g
102566	Papain, Technical powder. This is a crude, dried papaya latex which is not purified.	100 g 250 g 500 g 1 kg

References:

- Akahane, K., and Umeyama, H.: "Binding Specificity of Papain and Cathepsin B," *Enzyme*, **v. 36**, 141 (1986).
- Akalski, M., Lewis, S., Maggio, E., and Shafer, J.: "Effects of Proflavine on Papain's Catalytic Activity. A Proflavine-Mediated Decrease in K_m ," *Biochem.*, **v. 12**, 1884 (1973).
- Alecio, M., Dann, M., and Lowe, G.: "The Specificity of the S1 Subsite of Papain," *Biochem. J.*, **v.141**, 495 (1974).
- Allen, G., and Lowe, G.: "Investigation of the Active Site of Papain with Fluorescent Probes," *Biochem. J.*, **v. 133**, 679 (1973).
- Allison, W., and Swain, L.: "The Inactivation of Papain and Glyceraldehyde-3-Phosphate Dehydrogenase by Phenylidimide," *Arch. Biochem. Biophys.*, **v. 155**, 405 (1973).
- Angus, R., Carey, P., Lee, H., and Storer, A.: "Conformational Activation of Acylpapain and Acylcathepsin B's: Resonance Raman and Kinetic Evidence," *Biochem.*, **v. 25**, 3304 (1986).
- Arnon, R.: Papain, in *Methods in Enzymology*, **XIX**, (Perlmann, G., and Lorand, L., eds.), Academic Press, NY, 226 (1970).
- Arnon, R., and Shapira, E.: "Antibodies to Papain. A Selective Fractionation According to Inhibitory Capacity," *Biochem.*, **v. 6**, 3942 (1967).
- Asbóth, B., and Polgár, L.: "Transition-State Stabilization at the Oxyanion Binding Sites of Serine and Thiol Proteinases: Hydrolyses of Thiono and Oxygen Esters," *Biochem.*, **v. 22**, 117 (1983).
- Asbóth, B., Stokum E., Khan, I., and Polgár, L.: "Mechanism of Action of Cysteine Proteinases. Oxyanion Binding Site Is Not Essential in the Hydrolysis of Specific Substrates," *Biochem.*, **v. 24**, 606 (1985).
- Bader, C., MacLeish, P., and Schwartz, E.: "Responses to Light of Solitary Rod Receptors Isolated from Tiger Salamander

- Retina," *Proc. Natl. Acad. Sci. USA.*, **v. 75**, 3507 (1978).
- Balls, A., and Lineweaver, H.: "Isolation and Properties of Crystalline Papain," *J. Biol. Chem.*, **v. 130**, 669 (1939).
 - Banks, T., and Shafer, J.: "Inactivation of Papain by S-Methylation of Its Cysteiny Residue with O-Methylisourea," *Biochem.*, **v. 11**, 110 (1972).
 - Barrett, A., and Buttle, D.: "Names and Numbers of Papaya Proteinases," *Biochem. J.*, **v. 228**, 527 (1985).
 - Blumberg, S., Schechter, I., and Berger, A.: "The Purification of Papain by Affinity Chromatography," *Eur. J. Biochem.*, **v. 15**, 97 (1970).
 - Brisson, J., Carey, P., and Storer, A.: "Benzoylamidoacetonitrile Is Bound as a Thiomidate in the Active Site of Papain," *J. Biol. Chem.*, **v. 261**, 9087 (1986).
 - Brocklehurst, K., and Kierstan, M.: "Propapain and Its Conversion to Papain: A New Type of Zymogen Activation Mechanism Involving Intramolecular Thiol-Disulfide Interchange," *Nature New Biol.*, **v. 242**, 167 (1973).
 - Brocklehurst, K., and Little, G.: "Reactions of Papain and of Low-Molecular-Weight Thiols with Some Aromatic Disulfides," *Biochem. J.*, **v. 133**, 67 (1973).
 - Brocklehurst, K., and Salih, E.: "A Re-evaluation of the Nomenclature of the Cysteine Proteinases of Carica Papaya and a Rational Basis for their Identification," *Biochem. J.*, **v. 213**, 559 (1983).
 - Brocklehurst, K., Carey, P., Lee, H., Sahih, E., and Storer, A.: "Comparative Resonance Raman Spectroscopic and Kinetic Studies of Acyl-enzymes Involving Papain, Actinidin and Papaya Peptidase II," *Biochem. J.*, **v. 223**, 649 (1984).
 - Brocklehurst, K., Kierstan, M., and Little, G.: "The Reaction of Papain with Ellman's Reagent (5,5'-Dithiobis-2-nitrobenzoate Dianion)," *Biochem. J.*, **v. 128**, 811 (1972).
 - Brubacher, L., and Bender, M.: "The Preparation and Properties of Trans-Cinnamoyl-Papain," *J. Am. Chem. Soc.*, **v. 88**, 5871 (1966).
 - Brubacher, L., and Bender, M.: "The Reaction of Trans-Cinnamoyl-Papain with a Series of Polyglycinamides of Varying Chain Length," *Biochem. Biophys. Res. Comm.*, **v. 27**, 176 (1967).
 - Bryant, J., Leder, I., and Stetten, D.: "The Release of Chondroitin Sulfate from Rabbit Cartilage Following the Intravenous Injection of Crude Papain," *Arch. Biochem. Biophys.*, **v. 76**, 122 (1958).
 - Campbell, P., and Kaiser, E.: "The Reaction of A Cyclic Sulfonate Ester with the Sulfhydryl Proteolytic Enzyme Papain," *Bioorg. Chem.*, **v. 1**, 432 (1971).
 - Campbell, P., and Kaiser, E.: "The Intramolecular Nucleophilic Reaction of a Phenolic Hydroxyl Group at the Active Site of Papain," *Biochem. Biophys. Res. Comm.*, **v. 48**, 866 (1972).
 - Carotti, A., Smith, R., Wong, S., Hansch, C., Blaney, J., and Langridge, R.: "Papain Hydrolysis of X-Phenyl-N-methane-sulfonyl Glycinates: A Quantitative Structure-Activity Relationship and Molecular Graphics Analysis," *Arch. Biochem. Biophys.*, **v. 229**, 112 (1984).
 - Chaiken, I., and Smith, E.: "Reaction of a Specific Tyrosine Residue of Papain with Diisopropylfluorophosphate," *J. Biol. Chem.*, **v. 244**, 4247 (1969).
 - Clagett, J., and Englehard, W.: "Immunosuppressive Activity of the 3[[alpha]]-1 Fraction of Papain," *J. Bacteriol.*, **v. 93**, 1493 (1967).
 - Coates, J., and Swann, J.: "A Reinvestigation of the Association Behaviour of Papain in Solution," *Biochim. Biophys. Acta*, **v. 214**, 545 (1970).
 - Cohen, W., and Petra, P.: "A Kinetic Analysis of the Papain-Catalyzed Hydrolysis of [[alpha]]-N-Benzoyl-L-citrulline Methyl Ester," *Biochem.*, **v. 6**, 1047 (1967).
 - Coulson, W.: "The Effect of Proteolytic Enzymes on the Tensile Strength of Whole Aorta and Isolated Aortic Elastin," *Biochim. Biophys. Acta*, **v. 237**, 378 (1971).
 - deJersey, J.: "On the Specificity of Papain," *Biochem.*, **v. 9**, 1761 (1970).
 - Dreuth, J., Jansonius, J., Koekoek, R., Swen, H., and Wolters, B.: "Structure of Papain," *Nature*, **v. 218**, 929 (1968).
 - Fink, A., and Gwyn, C.: "Conformational Changes in Papain During Catalysis and Ligand Binding," *Biochem.*, **v. 13**, 1190 (1974).
 - Finkbeiner, S., and Stevens, C.F.: "Applications of Quantitative Measurements for Assessing Glutamate Neurotoxicity," *Proc. Nat. Acad. Sci. (USA)*, **v. 85**, 4071 (1988).
 - Frater, R., Light, A., and Smith, E.: "Chemical and Enzymic Studies on the Amino-Terminal Sequence of Papain. A Reinvestigation," *J. Biol. Chem.*, **v. 241**, 253 (1966).
 - Fujimoto, S., Kanazawa, H., Ishimitsu, S., and Ohara, A.: "On the Mechanism of Inactivation of Papain by Hydroxylamine," *Chem. Pharm. Bull. Tokyo*, **v. 38**, 546 (1990).
 - Glazer, A.: "Transesterification Reactions Catalyzed by Papain," *J. Biol. Chem.*, **v. 241**, 3811 (1966).
 - Glick, B., and Brubacher, L.: "Evidence for Nonproductive Binding Subsites Within the Active Site of Papain," *Can. J. Biochem.*, **v. 52**, 877 (1974).
 - Hall, P., Anderson, C., and Crawford, G.: "Acridine Dyes as Effectors in Papain Catalysis," *Arch. Biochem. Biophys.*, **v. 153**, 162 (1972).
 - Henry, A., and Kirsch, J.: "Papain-Catalyzed Reactions of Esters with Alcohols. The Nature of the Rate-Determining Step," *Biochem.*, **v. 6**, 3536 (1967).
 - Hill, R., and Smith, E.: *Biochim. Biophys. Acta*, **v. 19**, 376 (1956).
 - Hinkle, P., and Kirsch, J.: "Evidence for Conformational Restrictions within the Active Site of Acyl Papains Which Influence the Rates of Hydrolysis," *Biochem.*, **v. 9**, 4633 (1970).
 - Hinkle, P., and Kirsch, J.: "Evidence for Conformational Restrictions within the Active Site of Acyl Papains Which Influence the Rates of Hydrolysis. Correction," *Biochem.*, **v. 10**, 1510 (1971a).
 - Hinkle, P., and Kirsch, J.: "Kinetics and Thermodynamics of the Reactions of Acyl-Papains. Effect of pH, Temperature, Solvents, Ionic Strength and Added Nucleophiles," *Biochem.*, **v. 10**, 3700 (1971b).
 - Hu, L., and Abeles, R.: "Inhibition of Cathepsin B and Papain by Peptidyl [[alpha]]-Keto Esters, [[alpha]]-Keto Amides, [[alpha]]-Diketones, and [[alpha]]-Keto Acids," *Arch. Biochem. Biophys.*, **v. 281**, 271 (1990).
 - Huettner, J., and Baughman, R.: "Primary Culture of Identified Neurons from the Visual Cortex of Postnatal Rat," *J. Neuroscience*, **v. 6**, 3044 (1986).
 - Iso, N., and Williams, J.: "Sedimentation Analysis of a Papain-Digested Human Multiple Myeloma gamma-G-Globulin," *J. Biol. Chem.*, **v. 241**, 2787 (1966).
 - Jori, G., Gennari, G., Toniolo, C., and Scoffone, E.: "Probing the Topography of Proteins in Solution by Photosensitized Oxidation, The Catalytic Region of Papain," *J. Mol. Biol.*, **v. 59**, 151 (1971).
 - Kimmel, J., and Smith, E.: "Crystalline Papain. I. Preparation, Specificity and Activation," *J. Biol. Chem.*, **v. 207**, 515 (1954).

- Kimmel, J., Rogers, H., and Smith, E.: "Tryptic Digest of Papain. II. Isolation and Sequences of Peptides from the Carboxymethylated Protein," *J. Biol. Chem.*, **v. 240**, 266 (1965).
- Kirschenbaum, D.: "The Enhancement of the Enzymatic Activity of Papain by Reaction with N-Bromosuccinimide," *Biochim. Biophys. Acta*, **v. 235**, 159 (1971).
- Klein, I., and Kirsch, J.: "The Activation of Papain and the Inhibition of the Active Enzyme by Carbonyl Reagents," *J. Biol. Chem.*, **v. 244**, 5928 (1969).
- Kouvonen, I.: "Solubilization of the Pig Ileal Intrinsic Factor Receptor with Papain Treatment and Studies on the Solubilized Receptor," *Biochim. Biophys. Acta*, **v. 626**, 244 (1980).
- Kullmann, W.: "Kinetics of Chymotrypsin- and Papain-catalyzed Synthesis of [Leucine]enkephalin and [Methionine]enkephalin," *Biochem. J.*, **v. 220**, 405 (1984).
- Lake, A., and Lowe, G.: "The Kinetics of Papain- and Ficin-Catalyzed Hydrolyses in the Presence of Alcohols," *Biochem. J.*, **v. 101**, 402 (1966).
- Lam, D.: "Biosynthesis of Acetylcholine in Turtle Photoreceptors," *Proc. Natl. Acad. Sci. USA*, **v. 69**, 1987 (1972).
- Lewis, S., and Shafer, J.: "Ionization of a Nitrophenol-Containing Reporter Group at the Active Site of Papain," *Biochem.*, **v. 13**, 690 (1974).
- Liener, I.: *The Sulfhydryl Proteases, in Food Related Enzymes (Whitaker, J., ed.), Adv. in Chem. Series 136*, Amer. Chem. Soc., Washington DC, 202 (1974).
- Light, A., and Smith, E.: "Chymotryptic Digest of Papain. I.," *J. Biol. Chem.*, **v. 235**, 3144 (1960).
- Light, A., and Smith, E.: "Chymotryptic Digest of Papain. II.," *J. Biol. Chem.*, **v. 235**, 3151 (1960).
- Light, A., Frater, R., Kimmel, J., and Smith, E.: "Current Status of the Structure of Papain: The Linear Sequence, Active Sulfhydryl Group, and the Disulfide Bridges," *Proc. Natl. Acad. Sci.*, **v. 52**, 1276 (1964).
- Light, A., Glazer, A., and Smith, E.: "Chymotryptic Digest of Papain. III.," *J. Biol. Chem.*, **v. 235**, 3159 (1960).
- Lin, W., Armstrong, D., and Gaucher, G.: "Formation and Repair of Papain Sulfenic Acid," *Can. J. Biochem.*, **v. 53**, 298 (1975).
- Lindahl, P., Alriksson, E., Jörnvall, H., and Björk, I.: "Interaction of the Cysteine Proteinase Inhibitor Chicken Cystatin with Papain," *Biochem. USA*, **v. 27**, 5074 (1988).
- Lowe, G.: "Structural Relationships between Some Plant and Animal Proteases," *Nature*, **v. 212**, 1263 (1966).
- Lowe, G., and Yuthavong, Y.: "Kinetic Specificity in Papain-Catalyzed Hydrolyses," *Biochem. J.*, **v. 124**, 107 (1971a).
- Lowe, G., and Yuthavong, Y.: "pH-Dependence and Structure-Activity Relationships in the Papain-Catalyzed Hydrolysis of Anilides," *Biochem. J.*, **v. 124**, 117 (1971b).
- Luscombe, M., and Phelps, C.: "Action of Degradative Enzymes on the Light Fraction of Bovine Septa Protein Polysaccharide," *Biochem. J.*, **v. 103**, 103 (1967).
- Mackenzie, N., Grant, S., Scott, A., and Malthouse, J.: "¹³C NMR Study of the Sterospecificity of the Thiohemiacetals Formed on Inhibition of Papain by Specific Enantiomeric Aldehydes," *Biochem.*, **v. 25**, 2293 (1986).
- Martin, B., Feinman, R., and Detwiler, T.: "Plate Stimulation by Thrombin and Other Proteases," *Biochem.*, **v. 14**, 1308 (1975).
- McCluskey, R., and Thomas, L.: "The Removal of Cartilage Matrix, in vivo, by Papain. Identification of Crystalline Papain as the Cause of the Phenomenon," *J. Exp. Med.*, **v. 108**, 371 (1958).
- Migliorini, M., and Creighton, D.: "Active Site Ionizations of Papain - An Evaluation of the Potentiometric Difference Titration Method," *Eur. J. Biochem.*, **v. 156**, 189 (1986).
- Miller, W.: "Degradation of Poly-a, L-glutamic Acid. I. Degradation of High Molecular Weight PGA by Papain," *J. Am. Chem. Soc.*, **v. 83**, 259 (1961).
- Mitchel, R., Chaiden, I., and Smith, E.: "The Complete Amino Acid Sequence of Papain. Additions and Corrections," *J. Biol. Chem.*, **v. 245**, 3485 (1970).
- Mole, J., and Horton, H.: "2-Chloromethyl-4-Nitrophenyl(N-Carbobenzoxy)glucinate. A New Reagent Designed to Introduce an Environmentally Sensitive Conformational Probe Near the Active Site of Papain," *Biochem.*, **v. 12**, 5278 (1973a).
- Mole, J., and Horton, H.: "A Kinetic Analysis of the Enhanced Catalytic Efficiency of Papain Modified by 2-Hydroxy-5-Nitrobenzoylation," *Biochem.*, **v. 12**, 5285 (1973b).
- Mole, J., and Horton, H.: "Kinetics of Papain-Catalyzed Hydrolysis of a-N-Benzoyl-L-Arginine- p-Nitroanilide" *Biochem.*, **v. 12**, 816 (1973c).
- Morihara, K.: "Chemical Nature of the Active Site of Papain. I. Inhibition by Carbonyl Reagents with Active Methylene Groups," *J. Biochem. Japan*, **v. 62**, 250 (1967).
- Narinesingh, D., Ngo, T.: "Papain Covalently Immobilized on Fractogel Derivative: Preparation, Bioreactor Flow Kinetics, and Stability," *Biotech. Appl. Biochem.*, **v. 9**, 450 (1987).
- Neumann, H., Shinitzky, M., and Smith, R.: "The Activation of Papain and Ficin by Phosphorothioate," *Biochem.*, **v. 6**, 1421 (1967).
- Nicklin, M., and Barrett, A.: "Inhibition of Cysteine Proteinases and Dipeptidyl Peptidase I by Egg-white Cystatin," *Biochem. J.*, **v. 223**, 245 (1984).
- Nolan, C., and Smith, E.: "Absence of Sialic Acids and Other Carbohydrates in Crystalline Papain," *Proc. Soc. Exp. Biol. Med.*, **v. 105**, 287 (1960).
- Ozaki, Y., Pliura, D., Carey, P., and Storer, A.: "Vibrational Spectra of Scissile Bonds in Enzyme Active Sites: A Resonance Raman Study of Dithioacylpapains," *Biochem.*, **v. 21**, 3102 (1982).
- Ozawa, K., Ohnishi, T., and Tanaka, S.: "Activation and Inhibition of Papain," *J. Biochem.*, **v. 51**, 372 (1962).
- Polgár, L.: "Isolation of Highly Active Papaya Peptidases A and B from Commercial Chymopapain," *Biochim. Biophys. Acta*, **v. 658**, 262 (1981).
- Polgár, L.: "Problems of Classification of Papaya Latex Proteinases," *Biochem. J.*, **v. 221**, 555 (1984).
- Potter, J., McCluskey, R., Weissmann, G., and Thomas, L.: "The Removal of Cartilage Matrix by Papain," *J. Exp. Med.*, **v. 112**, 1173 (1960).
- Saunders, A., and Silverman, L.: "Electron Microscopy of Chondromucoprotein and the Products of its Digestion with Hyaluronidase and Papain," *Nature*, **v. 214**, 194 (1967).
- Schack, P., and Kaarsholm, N.: "Subsite Differences between Active Centres of Papaya Peptidase A and Papain as Revealed by Affinity Chromatography," *Biochem. J.*, **v. 219**, 727 (1984).
- Schack, P., and Kaarsholm, N.: "Absence in Papaya Peptidase A Catalyzed Hydrolysis of a pKa~4 Present in

- Papain-Catalyzed Hydrolyses," *Biochem.*, v. **23**, 631 (1984).
- Schechter, I., and Berger, A.: "On the Size of the Active Site in Proteases. I. Papain," *Biochem. Biophys. Res. Comm.*, v. **27**, 157 (1967).
 - Schneider, M., and Goldstein, G.: "Digestion of Articular Cartilage by Papain in vitro," *Proc. Soc. Exp. Biol. Med.*, v. **124**, 991 (1967).
 - Schwert, G., Neurath, H., Kaufman, S., and Snoke, J.: *J. Biol. Chem.*, v. **172**, 221 (1948).
 - Shapira, E., and Arnon, R.: "The Mechanism of Inhibition of Papain by Its Specific Antibodies," *Biochem.*, v. **6**, 3951 (1967a).
 - Shapira, E., and Arnon, R.: "Antibodies to Papain. Resolution on DEAE-Sephadex into Fractions with Different Precipitating Capacities," *Biochim. Biophys. Acta*, v. **140**, 177 (1967b).
 - Shipton, M., Kierstan, M., Malthouse, J., Stuchbury, T., and Brocklehurst, K.: "The Case for Assigning a Value of Approximately 4 to pKa1 of Essential Histidine-Cysteine Interactive Systems of Papain, Bromelain and Ficin," *FEBS Lett.*, v. **50**, 365 (1975).
 - Sluyterman, L.: "The Activation Reaction of Papain," *Biochim. Biophys. Acta*, v. **139**, 430 (1967).
 - Sluyterman, L.: "The Effect of Methanol, Urea, and Other Solutes on the Action of Papain," *Biochim. Biophys. Acta*, v. **139**, 418 (1967).
 - Sluyterman, L.: "Reversible Inactivation of Papain by Cyanate," *Biochim. Biophys. Acta*, v. **139**, 439 (1967).
 - Sluyterman, L., and DeGraaf, M.: "The Effect of Salts Upon the pH Dependence of the Activity of Papain and Succinyl-Papain," *Biochim. Biophys. Acta*, v. **258**, 554 (1972).
 - Sluyterman, L., and Wijdenes, J.: "Sigmoidal Progress Curves in the Polymerization of Leucine Methyl Ester Catalyzed by Papain," *Biochim. Biophys. Acta*, v. **289**, 194 (1972a).
 - Sluyterman, L., and Wijdenes, J.: "Cyanuration of Papain Activity and Fluorescence of the Products," *Biochim. Biophys. Acta*, v. **263**, 329 (1972b).
 - Sluyterman, L., and Wijdenes, J.: "Benzoylamido-acetonitrile as an Inhibitor of Papain," *Biochim. Biophys. Acta*, v. **302**, 95 (1973).
 - Smith, D., Maggio, E., and Kenyon, G.: "Simple Alkanethiol Groups for Temporary Blocking of Sulfhydryl Groups of Enzymes," *Biochem.*, v. **14**, 766 (1975).
 - Smith, E., and Kimmel, J.: in *The Enzymes*, Vol. **4**, 2nd Ed., (Boyer, P., Lardy, H., and Myrback, K., eds.), Academic Press, NY, 133 (1960).
 - Smith, E., Kimmel, J., and Brown, P.: *J. Biol. Chem.*, v. **207**, 533 (1954).
 - Smith, L., Garvin, P., Gesler, R., and Jennings, R.: "Enzyme Dissolution of the Nucleus Pulposus," *Nature*, v. **198**, 1311 (1963).
 - Steiner, R.: "Varying Luminescence Behavior of the Different Tryptophan Residues of Papain," *Biochem.*, v. **10**, 771 (1971).
 - Storer, A., and Carey, P.: "Comparison of the Kinetics and Mechanism of the Papain-Catalyzed Hydrolysis of Esters and Thiono Esters," *Biochem.*, v. **24**, 6808 (1985).
 - Storer, A., Lee, H., and Carey, P.: "Relaxed and Perturbed Substrate Conformations in Enzyme Active Sites: Evidence from Multichannel Resonance Raman Spectra," *Biochem.*, v. **22**, 4789 (1983).
 - Storer, A., Murphy, W., and Carey, P.: "The Use of Resonance Raman Spectroscopy to Monitor Catalytically Important Bonds during Enzymic Catalysis. Application to the Hydrolysis of Methyl Thionohippurate by Papain," *J. Biol. Chem.*, v. **254**, 3163 (1979).
 - Suzuki, H.: "Localized Hyperchondroplasia Induced by Crude Papain Administered to Mice," *Nature*, v. **209**, 1111 (1966).
 - Townes -Anderson, E., MacLeish, P., and Raviola, E.: "Rod Cells Dissociated from Mature Salamander Retina: Ultrastructure and Uptake of Horseradish Peroxidase," *J. Cell Biol.*, v. **100**, 175 (1985).
 - Ungar-Waron, H., Jatou, J., and Sela, M.: "Action of Papain on Normal Rabbit Immunoglobulin M," *Biochim. Biophys. Acta*, v. **140**, 542 (1967).
 - Valeri, A., Wilson, S., and Feinman, R.: "Reaction of Antithrombin with Proteinases. Evidence for a Specific Reaction with Papain," *Biochim. Biophys. Acta*, v. **614**, 526 (1980).
 - Van Eyk, H.: "Fragmentation of Human Globulin with Papain," *Biochim. Biophys. Acta*, v. **127**, 241 (1966).
 - Weinstein, B.: "The Evolution of Papain," *Biochem. Biophys. Res. Comm.*, v. **41**, 441 (1970).
 - Westerik, J., and Wolfenden, R.: "Aldehydes as Inhibitors of Papain," *J. Biol. Chem.*, v. **247**, 8195 (1972).
 - Whitehouse, M., and Leader, J.: "Biochemical Properties of Anti-Inflammatory Drugs. IX. Coupling of Oxidative Phosphorylation and Inhibition of a Thiol Enzyme (Papain) by Some Cyclic [[beta]]-Diones and Ninhydrin," *Biochem. Pharmacol.*, v. **16**, 537 (1967).
 - Williams, D., and Whitaker, J.: "Kinetics of Papain Catalyzed Hydrolyses of Neutral Substrates," *Biochem.*, v. **6**, 3711 (1967).
 - Wilson, A., and Gurner, B.: "Increased Affinity of Guinea Pig Thymocytes and Thymus-Dependent Lymphocytes for Papain-Treated Rabbit Erythrocytes Compared to Untreated Erythrocytes," *J. Immunol. Meth.*, v. **7**, 163 (1975).
 - Wolthers, B., and Kalk, K.: "Binding of Competitive (Product) Inhibitors to Papain in the Active and Inactive State," *Biochim. Biophys. Acta*, v. **198**, 556 (1970).
 - Wolthers, B., Dreuth, J., Jansonius, J., Koekoek, R., and Swen, H.: *The Three Dimensional Structure of Papain, in Structure-Function Relationships of Proteolytic Enzymes*, (Desuelle, P., Neurath, H., and Ottesen, M., eds.), Academic Press, NY, 272 (1970).
 - Zucker, S., Buttle, D., Nicklin, M., and Barrett, A.: "The Proteolytic Activities of Chymopapain, Papain and Papaya Proteinase III," *Biochim. Biophys. Acta*, v. **828**, 196 (1985).