



MP Biomedicals, LLC

29525 Fountain Parkway  
Solon, Ohio 44139

Telephone: 440/337-1200  
Toll Free: 800/854-0530  
Fax: 440/337-1180  
mailto: [biotech@mpbio.com](mailto:biotech@mpbio.com)  
web: <http://www.mpbio.com>

## TECHNICAL INFORMATION

Catalog Number: 55981  
**Whole Cat Serum**

**DESCRIPTION:** Product is whole cat serum and buffer salts.

The serum is delipidated, dialyzed into 0.02M sodium phosphate, 0.14M sodium chloride, pH 7.3, filtered through a 0.22 µm filter, vialled and lyophilized.

**USES:** This product is suitable for use as a serum protein control in immunoelectrophoresis, Ouchterlony immunodiffusion, and other immunoprecipitation techniques. It may also be used as a blocking agent in fluorescence and enzyme immunoassays.

**SPECIFICATIONS:** The total protein is measured using the Biuret procedure, with bovine albumin as the standard. The product gives a spectrum of normal serum precipitin arcs when tested by immunoelectrophoresis with anti-cat serum.

See **Related Product Ordering Information**.

**RECONSTITUTION:** Reconstitute product with 2.0 ml of deionized or distilled water. Gentle swirling may be used to speed rehydration. Avoid vigorous shaking of the reconstituted material.

### RELATED PRODUCT ORDERING INFORMATION:

Catalog #DescriptionSize

55288 (01031381)GOAT ANTI-CAT SERUM2ML  
55999 (51031381)GOAT ANTI-CAT SERUM (FORENSIC)2ML  
55830 (00001040)ID DISCS PATTERN-A WITH INDICATOR20EA  
55831 (00001050)ID DISCS PATTERN-A WITHOUT INDICATOR20EA

**Note:** This product may contain a preservative such as sodium azide, thimerosal or proclin. Please see lot specific chemical credential for preservative information.