



MP Biomedicals, LLC

29525 Fountain Parkway  
Solon, Ohio 44139

Telephone: 440/337-1200  
Toll Free: 800/854-0530  
Fax: 440/337-1180  
mailto: [biotech@mpbio.com](mailto:biotech@mpbio.com)  
web: <http://www.mpbio.com>

## TECHNICAL INFORMATION

# ImmunO™

Catalog Number: 623411

### **FITC-Streptavidin**

**Form:** Liquid in 0.01M sodium phosphate, 0.25M sodium chloride, pH 7.6, containing 15 mg/ml BSA, 40% glycerol. Streptavidin from *Streptomyces avidinii* is conjugated with FITC. The conjugate is purified by molecular sieving and is completely free of unbound FITC.

**Applications:** FITC-Streptavidin can be used for the detection or localization of biotinylated biomolecules.

**Working Concentration:** 2-5 ug/ml for most applications. Working dilution for specific applications should be determined by the investigator to obtain the best conditions.

**Note:** This product may contain a preservative such as sodium azide, thimerosal or proclin. Please see lot specific chemical credential for preservative information.