



MP Biomedicals, LLC

29525 Fountain Parkway
Solon, Ohio 44139

Telephone: 440/337-1200
Toll Free: 800/854-0530
Fax: 440/337-1180
mailto: biotech@mpbio.com
web: <http://www.mpbio.com>

TECHNICAL INFORMATION

ImmunO™

Catalog Number: 643491
Human IgE (Myeloma)

Description: Liquid. Product is in 15 mM potassium phosphate buffer, pH 7.4, with 0.15M NaCl.

Specificity: Single arc by IEP against antisera to whole human serum, human IgE and human kappa light chain.

Notes:

- By gel electrophoresis, it has been observed that when running the whole molecule of IgE, the heavy chain of myeloma IgE may appear as either a single or double band.
- Prepared from plasma shown to be nonreactive for HBsAg, anti-HCV, anti-HBc and negative for anti-HIV 1 and 2 by FDA-required tests.

Use Universal Precautions.

Note: This product may contain a preservative such as sodium azide, thimerosal or proclin. Please see lot specific chemical credential for preservative information.