



MP Biomedicals, LLC

29525 Fountain Parkway
Solon, Ohio 44139

Telephone: 440/337-1200
Toll Free: 800/854-0530
Fax: 440/337-1180
mailto: biotech@mpbio.com
web: <http://www.mpbio.com>

TECHNICAL INFORMATION

ImmunO™

Catalog Number: 650351
Rabbit Anti-Porcine Serum

Form: Lyophilized

Immunization Antigen: Normal pool of pig serum in Freund's adjuvant.

Immuno-electrophoresis: Versus porcine serum - multiple immunoprecipitin arcs observed after staining. The antiserum has not been assayed for interspecies cross reactivity.

Applications: Immuno-electrophoresis
Gel Diffusion Techniques

Note: This product may contain a preservative such as sodium azide, thimerosal or proclin. Please see lot specific chemical credential for preservative information.

[If a titer/working dilution is not given above, please click here to see a general dilution chart for working with antibodies. Please note that the general dilution chart should only be used as a guideline. Each lab should determine their own optimal working dilution.](#)

[Will this antibody work with your application? Please click here to see a general chart of antibody applications. Please note that any information given above is primary application data. The general applications charts should only be used as a reference.](#)