

Certificate of Analysis

Product: Agrobacterium Transformation Kit

Catalog No.: 3301-100

Lot No.: 103030

The Agrobacterium Transformation Kit contains a single transformation solution for use in either a freeze/thaw or electroporation method of plasmid DNA. Optimized Agrobacterium Medium and Agar are also included, along with Roll & Grow® Plating Beads for easy spreading of transformed cells. Once media preparation is complete, this easy transformation procedure takes less than 45 minutes with minimal hands-on time.

STORAGE CONDITIONS

All Agrobacterium Transformation Kits are shipped ambiently. The Agrobacterium Transformation Kit reagents are stable at room temperature. Expiration date is 05/31/2019.

SPECIFICATIONS

Reagent: Agrobacterium
Transformation Solution

Agrobacterium Transformation Solution

Cat #: 3301-101

Lot #: 101537

Passed

Reagent: Agrobacterium Media

Agrobacterium Media

Cat #: 3301-012

Lot #: 102688

Passed

Reagent: Agrobacterium Media Agar

Agrobacterium Media Agar

Cat #: 3301-212

Lot #: 102689

Passed

Reagent: Roll & Grow® Plating Beads

Roll & Grow® Plating Beads

Cat #: 5000-552

Lot #: 99649

Passed

REAGENT PERFORMANCE

All new lots of reagents are tested with the kit.

Conclusions:

This lot of Agrobacterium Transformation Kits are released for product sales.

MP Biomedicals certifies that this product was inspected and that all inspection operations and testing having been complete, the product conforms to the specifications.



07/05/2017 - John Huang, PhD
QBioGene/MP Biomedicals, LLC
Quality Control Manager

MSDS available online at www.mpbio.com