

Certificate of Analysis

Product Description: 4.0mm Glass Beads
Catalog Number: 6914YY0; 6914801; 2521011
Lot: 124048

Formula: N/A
CAS #: N/A
Physical Description: Clear, Colorless Glass Beads

Formula Weight: N/A
Storage: 15 - 30°C

Test	Specification	Result
Identity Test	Passes	Passes
DNase	Not Detected	Not Detected
RNase	Not Detected	Not Detected

Recommended Retest Date: 01/31/2029



01/08/2019 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control Manager

This is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)