

Certificate of Analysis

Product: GENECLEAN® DNA EasyClean MiniPrep Kit for Solid Tissues

Catalog No.: 6640-500; 6640500

Lot No.: 20396

The GENECLEAN® DNA EasyClean MiniPrep Kit for Solid Tissues is used for purification of DNA from solid samples.

STORAGE CONDITIONS


All GENECLEAN® DNA EasyClean MiniPrep Kits are shipped ambiently. GENECLEAN® DNA EasyClean MiniPrep Kit reagents are stable at room temperature.

SPECIFICATIONS

Component	Catalog #	Lot #
Proteinase K, 20 mg	6640302	20392
Proteinase K Mixx, 1.2 ml	6640202	20387
Digestion Buffer, 20 ml	6610301	20391
GTC Lysis Buffer, 100 ml, UN1760-III, Class 8	6580301	20373
IPA Pre-Wash Buffer, 50 ml, UN1219-II, Class 3	6580305	20372
DNA Wash Buffer, 100 ml, UN1219-II, Class 3	6580302	20371
DNA Elution Buffer, 50 ml	6610304	20390

MP Biomedicals, LLC certifies that this product was inspected and that all inspection operations and testing having been complete, the product conforms to the specifications.

Recommended Retest Date: 6/11/2017



06/12/2012 - Christina Marotta
QBioGene/MP Biomedicals, LLC
Technical Director

MSDS available online at www.mpbio.com