



MP Biomedicals, LLC

29525 Fountain Parkway

Solon OH 44139

United States

800.854.0530

www.mpbio.com

Certificate of Analysis

Description: Cell Culture Flask, 25 cm ² , vent cap, TC-treated	
Catalog number: 096140200	QI Request No: 00006961
Related Catalog number: 096140200R	
Lot number: 2301001C	Quality Level: N/A
Recommended Retest Date (MM/DD/YYYY): 01-24-2028	

Formula: N/A	Formula Weight: N/A
CAS #: N/A	Storage: 15-30 °C
	Storage Conditions: N/A

Additional Product information	Additional Lot Information
LOD: Limit of detection	

Test Description	Specification	Result	Comments
Identity Test	PASS	PASS	
DNase detection (LOD: 1 x10 ⁻⁷ Kunitz units) - Not Detected	PASS	PASS	
Physical Description (Cell Culture Flask)	PASS	PASS	Cell Culture Flask
RNase Detection (LOD: 1 x10 ⁻⁹ Kunitz units) - Not Detected	PASS	PASS	

05-16-2023 Fengqi An
MP Biomedicals, LLC.

This is an electronically generated document

Mail to: biotech@mpbio.com