

Certificate of Quality

1.0 Product Name: Smart Vortenizer

2.0 Catalogue no: 118100500

3.0 Lot no: S2604001010

4.0 Serial no: 2604018C

5.0 Parameters Checked:

5.1 Visual Inspection:

5.1.1 The unit and packaging box are free from damage.

5.1.2 The serial number on unit and packaging box are identical.

5.1.3 The red light illuminates when powered on.

5.2 Functional Test:

5.2.1 Speed Increment & Decrement Function Test: Pass

5.2.2 Time Increment & Decrement Function Test : Pass

5.2.3 LED Function Test (Progress bar) : Pass

5.2.4 Countdown and Stopwatch Function Test : Pass

5.2.5 Sensor Mode Function Test : Pass

5.2.6 Machine Noise Level Measurement : Pass

The above inspection criteria's have been reviewed and approved according to the Quality Assurance Procedures of MP Biomedical (Shandong) Co. Ltd.

Chris. Qu

Chris Qu
Quality manager

Date:

2026-07-28

