

Certificate of Analysis

Product: Salton Wash #2
Catalog No.: 1002-208
Lot No.: 29780

SPECIFICATIONS

Reagent:	Proprietary	Salton	Wash	#2									
pH	Range:	7.0-7.4	Cat	#:	1002-208								
6	M	solution	Lot		#:29780								
Reagent	melted	0.1 g	gel	slice	in	5	minutes	at	55°C	Passed			
Less	than	15%	DNA	1	kb	and	above	is	left	unbound	in	the	reagent.


REAGENT

Overexposure to oxygen or contaminants will cause oxidation and subsequently a rise in the pH of the solution. Thus reducing the shelf life. Yellowing of solution may be an indication of such oxidation.

PERFORMANCE

Conclusions:

This lot of Salton Wash #2 is released for product sales.



01/03/2013 - Christina Marotta
QBioGene/MP Biomedicals, LLC
Technical Director

MSDS available online at www.mpbio.com