



MP Biomedicals, LLC

29525 Fountain Parkway

Solon OH 44139

United States

800.854.0530

www.mpbio.com

Certificate of Analysis

Description: SuperFastPrep-2™	
Catalog number: 116012500	QI Request No: 00037280
Related Catalog number:	
Lot number: 427	Quality Level: N/A
Recommended Retest Date (MM/DD/YYYY): 10-31-2035	

Formula: N/A	Formula Weight: N/A
CAS #: N/A	Storage: 15 - 30°C
	Storage Conditions: N/A

Additional Product information	Additional Lot Information
---------------------------------------	-----------------------------------

Test Description	Specification	Result	Comments
Identity	PASS	PASS	
Capacity (2 x 2 mL Tube)	PASS	PASS	
Dimensions (13"L x 3.4"W x 4.6"H)	PASS	PASS	
Physical Description (SuperFastPrep-2)	PASS	PASS	
Sample Volume (up to 200 mg in 1 ml buffer)	PASS	PASS	

04-27-2026 Blessing Enya
MP Biomedicals, LLC.

This is an electronically generated document

Mail to: biotech@mpbio.com