

Certificate of Analysis

Product: FastPROTEIN RED (50 prep)
Catalog No.: 6550-600
Lot No.: 57124

The FastPROTEIN Red Kit is designed for rapid isolation and purification of protein from yeast cultures. The process is designed to be a small scale test for protein induction.

STORAGE CONDITIONS:

All FastPROTEIN Red Kits are shipped at ambient temperature. Expiration Date: August, 2017. Store YBB at 4°C.

SPECIFICATIONS:

Reagent:

YBB
pH Range: 7.2 - 7.6
Made the cells homogenous without lysing the cells.

YBB
Cat #: 6550-601
Lot #: 56709
Passed

Lysing Matrix C


Lysing Matrix C
Cat #: 6912-050
Lot #: 51890
Passed

REAGENT PERFORMANCE:

All reagents are stable at room temperature

In order to pass quality control tests:

The yield of the protein obtained using the new reagents must be at least 10ug/1.5ml sample and within 10% of the yield of the test standard lot.



08/05/2014 - Christina Marotta
QBioGene/MP Biomedicals, LLC
Technical Director

MSDS available online at www.mpbio.com