

Certificate of Analysis

Product: Agrobacterium Transformation Kit
Catalog No.: 3301-100
Lot No.: 68836

The Agrobacterium Transformation Kit contains a single transformation solution for use in either a freeze/thaw or electroporation method of plasmid DNA. Optimized Agrobacterium Medium and Agar are also included, along with Roll & Grow® Plating Beads for easy spreading of transformed cells. Once media preparation is complete, this easy transformation procedure takes less than 45 minutes with minimal hands-on time.

STORAGE CONDITIONS

All Agrobacterium Transformation Kits are shipped ambiently. The Agrobacterium Transformation Kit reagents are stable at room temperature. Expiration date is 03/31/2018.

SPECIFICATIONS

Reagent: Agrobacterium
Transformation Solution

Agrobacterium Transformation Solution
Cat #: 3301-101
Lot #: 68782
Passed

Reagent: Agrobacterium Media

Agrobacterium Media
Cat #: 3301-025
Lot #: 67983
Passed

Reagent: Agrobacterium Media Agar

Agrobacterium Media Agar
Cat #: 3301-225
Lot #: 67984
Passed

Reagent: Roll & Grow® Plating Beads

Roll & Grow® Plating Beads
Cat #: 5000-552
Lot #: 68858
Passed

REAGENT PERFORMANCE

All new lots of reagents are tested with the kit.

Conclusions:

This lot of Agrobacterium Transformation Kits are released for product sales.

MP Biomedicals certifies that this product was inspected and that all inspection operations and testing having been complete, the product conforms to the specifications.



06/17/2015 - John Huang, PhD
QBioGene/MP Biomedicals, LLC
Senior QC Chemist

MSDS available online at www.mpbio.com