

Certificate of Analysis

Product: FastPROTEIN RED (100 prep)
Catalog No.: 6550-700
Lot No.: 79318

The FastPROTEIN Rd Kit is designed for rapid isolation and purification of protein from yeast cultures. The process is designed to be a small scale test for protein induction.

STORAGE CONDITIONS:

All FastPROTEIN Red Kits are shipped at ambient temperature. Expiration Date: November, 2018. Store YBB at 4°C.

SPECIFICATIONS:

Reagent: YBB
ph Range: 7.2-7.6
Made the cells homogenous without lysing the cells.

YBB
Cat #: 6550-701
Lot #: 78987
Passed

Lysing Matrix C

Lysing Matrix C
Cat #: 6912-100
Lot #: 79146
Passed

REAGENT PERFORMANCE:

All reagents are stable at room temperature

In order to pass quality control tests:

The yield of the protein obtained using the new reagents must be at least 10ug/1.5ml sample and within 10% of the yield of the test standard lot.



12/07/2015 - John Huang, PhD
QBioGene/MP Biomedicals, LLC
Quality Control Manager

MSDS available online at www.mpbio.com