

Certificate of Analysis

Product: CLS-TC
Catalog No.: 6540-X09; 6540-409
Lot No.: 84296

STORAGE CONDITIONS

This reagent is stable at room temperature. Shelf life can be extended by storing at 4°C.

SPECIFICATIONS

CLS-TC is a colorless, precipitate-free solution, and has a pH range of 5.8 - 6.2.

Test	Specification	Result
Clarity	Colorless, precipitate-free solution	Passes
pH	5.8 - 6.2	6.11
Expiration Date	Month/Day/Year	03/31/2020



03/30/2016 - John Huang, PhD
QBioGene/MP Biomedicals, LLC
Quality Control Manager

MSDS available online at www.mpbio.com