



MP Biomedicals, LLC
29525 Fountain Parkway

Solon OH 44139
United States
800.854.0530
www.mpbio.com

Certificate of Analysis

Description: MagBeads Fast Circulating DNA Kit (Ready-to-Use for MPure-32™)

Catalog number: 117033900

QI Request No: N/A

Related Catalog number: N/A

Lot number: FE01035

Quality Level: N/A

Expiration Date (MM/DD/YYYY): 06-08-2025

Formula: N/A

Formula Weight: N/A

CAS #: N/A

Storage: 15 - 30°C

Storage Conditions: N/A

Additional Product information

Additional Lot Information

Test Description	Specification	Result	Comments
Appearance of the solutions	Complete components, clean appearance, no leakage, no damage. Colorless and transparent solution. The Magnetic Beads is black suspension.	PASS	
Physical and chemical properties of solution	The density, conductivity and pH of the solution meet the requirements	PASS	

Julie Ma

01-21-2024 *Julie Ma*
MP Biomedicals, LLC.

This is an electronically generated document

Mail to: biotech@mpbio.com