

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Certificate of Analysis

Product Description: Poly-L-Lysine Hydrobromide, cell culture reagent, gamma -irradiated
Catalog Number: 194543
Lot: MR29016

Formula: N/A
CAS #: 25988-63-0
Physical Description: Off-white Flaky Powder

Formula Weight: 30,000 - 70,000 kDa
Storage: -0°C

Test	Specification	Result
Identity Test	Passes	Passes
MW	30,000 - 70,000 kD	30,000-65,000 Da (by viscosity) 30,000 - 70,000 Da (by HPLC)
gamma-Irradiated	Passes	Passes

Source: Synthetic
Counter ion: HBr
Purity (NMR): 95%

03/27/2012 - Christina Marotta
MP Biomedicals, LLC.
Technical Director

This is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)