

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Certificate of Analysis

Product Description: Phosphatidyl-L-serine
Catalog Number: 100969
Lot: Q3597**Formula:** N/A
CAS #: N/A
Physical Description: Clear, Pale Yellow Solution**Formula Weight:** 788
Storage: -20°C

Test	Specification	Result
Identity Test	Passes	Passes
Purity	~98%	>98% (TLC)
Source	Bovine Brain	Bovine Brain
Concentration	10 - 40 mg/ml	10 mg/ml
Supplied in	Chloroform	Chloroform

Remove solvent under vacuum.

Fatty acid: Primarily stearic acid (C18:0) and Oleic (C18:1) acids

Recommended Retest Date: 07/06/2017

02/04/2016 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control ManagerThis is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)