

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Certificate of Analysis

Product Description: Prostaglandin e2
Catalog Number: 195461
Lot: Q4365**Formula:** $C_{20}H_{32}O_5$
CAS #: 363-24-6
Physical Description: White Solid**Formula Weight:** 352.5
Storage: -20°C, desiccate

Test	Specification	Result
Identity Test	Passes	Passes
Purity	≥99%	100% (HPLC)
Source	Synthetic	Synthetic

Mass/Source type:

- M-H⁺: 351.9
- H₂O: 333.4
- -2H₂O: 315.5
- -2H₂O, CO₂: 271.6
- CH₃ (CH₂)₄CHO, H₂O: 233.0

Melting point: 68 - 69°C

Purity (TLC): 100.00%

Rf: 0.44

Recommended Retest Date: 02/09/2018

A handwritten signature in black ink, appearing to read 'J Huang'.

06/15/2016 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control ManagerThis is an electronically generated document
mailto:biotech@mpbio.com
http://www.mpbio.comOnline Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)