



Certificate of Analysis

Product Description: 15 ml Conical Tubes (500)
Catalog Number: DP1090B
Lot: Q4503

Formula: N/A
CAS #: N/A
Physical Description: Clear Plastic Conical Tube, with white space for writing, Blue Cap

Formula Weight: N/A
Storage: Room Temperature

Test	Specification	Result
Identity Test	Passes	Passes

Recommended Retest Date: 7/10/2021

A handwritten signature in black ink, appearing to read 'J Huang'.

07/11/2016 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control Manager

This is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)