

**MP Biomedicals, LLC**29525 Fountain Parkway  
Solon, Ohio 44139Telephone: 440/337-1200  
Toll Free: 800/854-0530

## Certificate of Analysis

**Product Description:** Rabbit, Anti-Bovine Albumin  
**Catalog Number:** 65111  
**Lot:** Q8365

**Formula:** N/A  
**CAS #:** N/A  
**Physical Description:** Brown Liquid

**Formula Weight:** N/A  
**Storage:** 2 - 8°C

Test	Specification	Result
Identity Test	Passes	Passes
Host	Rabbit	Rabbit
Form	Lyophilized	Liquid, Need to Lyophilize
Antigenic Specificity	Bovine Albumin	Bovine Albumin
Concentration	For Information Only	56 mg/ml

Recommended Retest Date: 02/07/2019



01/16/2018 - John Huang, PhD  
MP Biomedicals, LLC.  
Quality Control Manager

This is an electronically generated document  
<mailto:biotech@mpbio.com>  
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site  
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)