

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530Fax: 440/337-1180
web: www.mpbio.com

Certificate of Analysis

Product Description: Albumin, Anti-Bovine, (RB) IgG
Catalog Number: 55285
Lot: Q8485**Formula:** N/A
CAS #: N/A
Physical Description: Peach to Brown Solid**Formula Weight:** N/A
Storage: -20°C

Test	Specification	Result
Identity Test	Passes	Passes
Purity	≥ 95%	≥ 95% (SDS-Page)
Concentration	Report Result	30 mg/ml
Formulation	PBS (no preservative)	PBS (no preservative)

Specificity:

- Bovine Albumin Purified Protein: one band observed
- Normal Bovine Serum: one band observed
- Bovine IgG and Bovine Transferrin: negative for both

ELISA: Expected range of 1:32K - 1:64K
Recommended Retest Date: 01/30/202101/31/2018 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control ManagerThis is an electronically generated document
mailto:biotech@mpbio.com
http://www.mpbio.comOnline Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)