

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Certificate of Analysis

Product Description: I-HUMAN FIBRONECTIN
Catalog Number: I991010830
Lot: Q8607**Formula:** N/A
CAS #: N/A
Physical Description: White Lyophilized Solid**Formula Weight:** 550000
Storage: 2 - 8°C

Test	Specification	Result
Identity Test	Passes	Passes
Concentration	Report Result	Lyophilized@1.7 mg/ml
Recommended Retest Date	Month/Day/Year	02/15/2019

Buffer: 20 mM Hepes, 150 mM NaCl, pH 7.4
E/1% (1cm, 280nm) = 14.002/20/2018 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control ManagerThis is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)