

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Certificate of Analysis

Product Description: I-PA-RB-A-BOV-ALBUMIN
Catalog Number: I996020342
Lot: Q8681

Formula: N/A
CAS #: N/A
Physical Description: Clear, Colorless Liquid

Formula Weight: N/A
Storage: 2 - 8°C

| Test | Specification | Result |
|---------------|-----------------------------------|-----------------------------------|
| Identity Test | Passes | Passes |
| Host | Rabbit | Rabbit |
| Iso Type | IgG | IgG |
| Buffer | Phosphate buffered saline, pH 7.2 | Phosphate buffered saline, pH 7.2 |

Preservative: 0.1% sodium azide
Purity: Antigen Affinity Purified
Recommended Retest Date: 02/28/2020



03/05/2018 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control Manager

This is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)