

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Certificate of Analysis

Product Description: Linolenic Acid, 45-65%
Catalog Number: 102188
Lot: Q8924**Formula:** N/A
CAS #: N/A
Physical Description: Clear, Yellow Liquid**Formula Weight:** N/A
Storage: -20°C

Test	Specification	Result
Identity Test	Passes	Passes
Purity	45 - 65%	53.6%
FTIR	Conforms to Reference	Conforms to Reference

A natural triglyceride oil containing about 50% linolenic acid.

Peroxide Value: 2.9 meq/kg
Heavy metals: None detected

Fatty acid profile:

- Palmitic: 5.09 area%
- Stearic: 3.8 area%
- Oleic: 22 area%
- Linoleic: 14.6 area%

Recommended Retest Date: 08/31/2020



04/04/2018 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control Manager

This is an electronically generated document
mailto:biotech@mpbio.com
http://www.mpbio.com

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)