

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Certificate of Analysis

Product Description: Tributyrin
Catalog Number: 103111
Lot: Q9351**Formula:** $C_{15}H_{26}O_6$
CAS #: 60-01-5
Physical Description: Clear, Colorless Liquid**Formula Weight:** 302.36
Storage: 15 - 30°C

Test	Specification	Result
Identity Test	Passes	Passes
Purity	≥96%	96.9%
FTIR	Conforms to Standard	Conforms to Standard

Water: ≤0.1%
Acid Value: ≤5.0
Refractive Index @ 20°C: 1.435
Specific Gravity @ 25°C: 1.031
Solubility: Passes Test
Recommended Retest Date: 05/15/2020



06/05/2018 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control Manager

This is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)