

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Certificate of Analysis

Product Description: MP Oakton® Conductivity/TDS Solution, 23 µS
Catalog Number: 160145
Lot: QR12760**Formula:** N/A
CAS #: N/A
Physical Description: Clear, Colorless Liquid**Formula Weight:** N/A
Storage: Room Temperature

Test	Specification	Result
Identity Test	Passes	Passes

Conductivity (25°C, traceable to NIST SRM 3193): 23.1 µS/cm
Recommended Retest Date: 08/15/201809/16/2016 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control ManagerThis is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)