

Certificate of Analysis

Product Description: ARIES Stool Resuspension Tubes
Catalog Number: 6911030
Lot: QR12909

Formula: N/A
CAS #: N/A
Physical Description: White Glass Beads (Silica)

Formula Weight: N/A
Storage: Room Temperature

Test	Specification	Result
Identity Test	Passes	Passes
Size Range	79-106 microns	Passes
Sphericity	≥ 90%	Passes
Heating Conditions	for 12-18 hours at 240-250°C	Passes
DNase	Not Detected	Not Detected
RNase	Not Detected	Not Detected

Specific gravity 2.46-2.49 g/cm³
Recommended Retest Date: 9/30/2026



10/14/2016 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control Manager

This is an electronically generated document
mailto:biotech@mpbio.com
http://www.mpbio.com

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)