

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Certificate of Analysis

Product Description: Succinic Anhydride
Catalog Number: 102970
Lot: QR14090**Formula:** C₄H₄O₃
CAS #: 108-30-5
Physical Description: White flakes**Formula Weight:** 100.1
Storage: 15-30°C desiccate

Test	Specification	Result
Identity Test	Passes	Passes
Purity	≥99%	99.6%
FTIR	Conforms to Standard	Conforms to Standard

Color, Molten (Pt-Co): 120
Recommended Retest Date: 05/10/202105/17/2017 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control ManagerThis is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)