

Certificate of Analysis

Product Description: GeneClean DNA EasyClean MiniPrep Kit solid tissues
Catalog Number: 6640500; 6640-500
Lot: QR14755

Formula: N/A
CAS #: N/A
Physical Description: Kit

Formula Weight: N/A
Storage: 15 - 30°C

The GENE CLEAN® DNA EasyClean MiniPrep Kit for Solid Tissues is used for purification of DNA from solid samples.

Component	Catalog #	Lot #
Proteinase K, 20 mg	6640302	20392
Proteinase K Mixx, 1.2 ml	6640202	20387
Digestion Buffer, 20 ml	6610301	20391
GTC Lysis Buffer, 100 ml, UN1760-III, Class 8	6580301	20373
IPA Pre-Wash Buffer, 50 ml, UN1219-II, Class 3	6580305	20372
DNA Wash Buffer, 100 ml, UN1219-II, Class 3	6580302	20371
DNA Elution Buffer, 50 ml	6610304	20390

Test	Specification	Result
Kit Test	Passes	Passes

MP Biomedicals, LLC certifies that this product was inspected and that all inspection operations and testing having been complete, the product conforms to the specifications.

Recommended Retest Date: 09/19/2019



10/12/2017 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control Manager

This is an electronically generated document
mailto:biotech@mpbio.com
http://www.mpbio.com

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)