

**MP Biomedicals, LLC**29525 Fountain Parkway  
Solon, Ohio 44139Telephone: 440/337-1200  
Toll Free: 800/854-0530

## Certificate of Analysis

**Product Description:** Goat Anti-Rabbit IgG (Whole Molecule), Affinity Pure, Alkaline Phosphatase Conjugated  
**Catalog Number:** 59298  
**Lot:** S0603

**Formula:** N/A  
**CAS #:** N/A  
**Physical Description:** Solid

**Formula Weight:** N/A  
**Storage:** 2 - 8°C

Test	Specification	Result
Identity Test	Passes	Passes
Host	Goat	Goat
Form	Lyophilized powder, Affinity Pure	Lyophilized powder, Affinity Pure
Antigenic Specificity	Rabbit IgG (whole molecule)	Rabbit IgG (whole molecule)
Conjugate	Alkaline Phosphatase	Alkaline Phosphatase

Buffer: 0.01M Tris-HCl, 0.25M NaCl, pH 8.0

Stabilizer: 15mg/mL Bovine Serum Albumin (IgG-Free, Protease-Free)

Preservative: 0.05% Sodium Azide

Suggested dilution range: 1:5,000 - 1:50,000 for ELISA and Western blotting

Recommended Retest Date: 12/19/2023

12/20/2018 - John Huang, PhD  
MP Biomedicals, LLC.  
Quality Control Manager

This is an electronically generated document  
mailto:biotech@mpbio.com  
http://www.mpbio.com

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site  
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)