

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Certificate of Analysis

Product Description: Anti-Human IgG Monoclonal Antibody from Mouse
Catalog Number: 63481
Lot: S1087

Formula: N/A
CAS #: N/A
Physical Description: Clear, Colorless Liquid

Formula Weight: N/A
Storage: 2 - 8°C

Test	Specification	Result
Identity Test	Passes	Passes
Concentration	Report Result	7.1 mg/ml

Buffer: 10mM Phosphate Buffer (at pH7.4) containing 0.9% NaCl and 0.02% NaN3.
Avoid repeated freezing and thawing
Recommended Retest Date: 05/31/2022



03/06/2019 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control Manager

This is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)