

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Certificate of Analysis

Product Description: Rabbit, Anti-Hemagglutinin (HA), Affinity Pure
Catalog Number: 55138
Lot: S2039

Formula: N/A
CAS #: N/A
Physical Description: Clear Colorless Liquid

Formula Weight: N/A
Storage: 2 - 8°C

Test	Specification	Result
Identity Test	Passes	Passes
Host	Rabbit	Rabbit
Form	Liquid, Affinity Pure	Liquid, Affinity Pure
Immunogen	Influenza Hemagglutinin-HA-Epitope (YPYDVPDYA) Conjugated to KLH	Influenza Hemagglutinin-HA-Epitope (YPYDVPDYA) Conjugated to KLH
Antigenic Specificity	Hemagglutinin	Hemagglutinin

Recommended Retest Date (MM/DD/YYYY): 07/01/2021

07/03/2019 - John Huang, PhD
 MP Biomedicals, LLC.
 Quality Control Manager

This is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
 Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)