

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Certificate of Analysis

Product Description: I-IgG-GT-A-Cat-IgG (H&L)
Catalog Number: I992030081
Lot: S2184

Formula: N/A
CAS #: N/A
Physical Description: Frozen Liquid

Formula Weight: N/A
Storage: -20°C

Test	Specification	Result
Identity Test	Passes	Passes
Protein Concentration	Report Result	25mg/mL
Titer	Report Result	≥ 1.15mg/ml
Specificity	Cat IgG Heavy and Light Chain	Cat IgG Heavy and Light Chain
Immunogen	Highly Purified Cat IgG	Highly Purified Cat IgG

Preservative: 0.1% NaN₃

Avoid repeated freeze/thaws

Recommended Retest Date (MM/DD/YYYY): 07/01/2020



07/26/2019 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control Manager

This is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)