

Certificate of Analysis

Product Description: FastGlycoProtein Isolation Kit WGA Resin
Catalog Number: 6550900
Lot: S2185

Formula: N/A
CAS #: N/A
Physical Description: 10 Preps Kit

Formula Weight: N/A
Storage: 2 - 8°C

The FastGlycoProtein™ Isolation Kit WGA Resin can be used for isolating proteins with sialic acid/GlcNAC-containing glycans.

Component	Cat#	Lot#	UN#, Packing Group, Classification	Volume
WGA Lectin Resin	6550804	TE269345		2.2 mL
Binding/Wash Buffer, 5X	6550802	TG270501		6.5 mL
Elution Buffer	6550803	TE266794		5 mL

Test	Specification	Result
Identity Test	Passes	Passes

The FastGlycoProtein™ Isolation Kit WGA Resin is shipped on cold packs. The product should be stored refrigerated upon arrival.

Recommended Retest Date (MM/DD/YYYY): 07/24/2023



07/25/2019 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control Manager

This is an electronically generated document
mailto:biotech@mpbio.com
http://www.mpbio.com

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)