

**MP Biomedicals, LLC**29525 Fountain Parkway  
Solon, Ohio 44139Telephone: 440/337-1200  
Toll Free: 800/854-0530

## Certificate of Analysis

**Product Description:** Mouse, Anti-Epithelial Keratin, AE 5  
**Catalog Number:** 69143  
**Lot:** S2671

**Formula:** N/A  
**CAS #:** N/A  
**Physical Description:** Frozen Clear Liquid

**Formula Weight:** N/A  
**Storage:** -20°C

Test	Specification	Result
Identity Test	Passes	Passes
Source	Mouse	Mouse
Form	Liquid, IgG Fraction	Liquid, IgG Fraction
Clone	AE5	AE5
Antigenic Specificity	Epithelial Keratin, 64 kDa	Epithelial Keratin, 64 kDa
Concentration	For Information Only	1mg/ml

Recommended Retest Date (MM/DD/YYYY): 10/01/2024



10/03/2019 - John Huang, PhD  
MP Biomedicals, LLC.  
Quality Control Manager

This is an electronically generated document  
<mailto:biotech@mpbio.com>  
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site  
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)