

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Certificate of Analysis

Product Description: Peroxidase-conjugated Goat Affinity Purified Antibody to Rabbit IgG (Whole Molecule)
Catalog Number: 55689
Lot: S3089

Formula: N/A
CAS #: N/A
Physical Description: Solid

Formula Weight: N/A
Storage: 2 - 8°C

Test	Specification	Result
Identity Test	Passes	Passes
Antibody Concentration	Report Result	0.8 mg/ml
Buffer	Report Result	0.01M Sodium Phosphate
Stabilizer	Report Result	15 mg/ml Bovine Serum Albumin (IgG-Free, Protease Free)
Antibody Specificity	Reacts with Whole Molecule of Rabbit IgG	Reacts with Whole Molecule of Rabbit IgG

Preservative: None

Peroxidase: Purified from horseradish root

Purity: The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads

Reconstitution: When ready to use, rehydrate with 2.0ml DI H₂O and centrifuge if not clear.

Recommended Retest Date (MM/DD/YYYY): 11/26/2024

06/04/2020 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control Manager

This is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)