

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Fax: 440/337-1180

web: www.mpbio.com

Certificate of Analysis

Product Description: Goat Anti-Rabbit IgG (Whole Molecule), Affinity Pure, Alkaline Phosphatase Conjugated
Catalog Number: 59298
Lot: S3206

Formula: N/A
CAS #: N/A
Physical Description: Powder

Formula Weight: N/A
Storage: 2 - 8°C

Test	Specification	Result
Identity Test	Passes	Passes
Host	Goat	Goat
Form	Lyophilized powder, Affinity Pure	Lyophilized powder, Affinity Pure
Antigenic Specificity	Rabbit IgG (whole molecule)	Rabbit IgG (whole molecule)
Conjugate	Alkaline Phosphatase	Alkaline Phosphatase

Buffer: 0.01M Tris-HCl, 0.25M NaCl, pH 8.0

Stabilizer: 15mg/mL Bovine Serum Albumin (IgG-Free, Protease-Free)

Preservative: 0.05% Sodium Azide

Suggested dilution range: 1:5,000 - 1:50,000 for ELISA and Western blotting

Recommended retest date (MM/DD/YYYY): 12/15/2024

A handwritten signature in black ink that reads 'John Huang'.

12/17/2019 - John Huang, PhD

MP Biomedicals, LLC.

Quality Control Manager

This is an electronically generated document

<mailto:biotech@mpbio.com><http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site

Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)